

Mock-community	Name	Aerobe/anaerobe	Temperature (°C)	Cultivation media
ZIEL-I-mock	<i>Actinomyces bowdenii</i>	aerobe	37	TSA/TSB
	<i>Enterorhabdus mucosicola</i>	anaerobe	37	WCA
	<i>Cellulosimicrobium cellulans</i>	anaerobe	37	WCA
	<i>Bacteroides sartorii</i>	anaerobe	37	WCA
	<i>Alistipes</i> sp.	anaerobe	37	WCA
	<i>Bacillus subtilis</i>	aerobe	37	TSA/TSB
	<i>Parabacteroides goldsteinii</i>	anaerobe	37	WCA
	<i>Flavonifractor plautii</i>	anaerobe	37	WCA
	<i>Clostridium ramosum</i>	anaerobe	37	WCA
	<i>Enterococcus hirae</i>	aerobe	37	TSA/TSB
	<i>Acetatifactor muris</i>	anaerobe	37	WCA
	<i>Staphylococcus warneri</i>	aerobe	37	TSA/TSB
	<i>Pseudomonas</i> sp.	aerobe	37	TSA/TSB
ZIEL-II-mock	<i>Prevotella copri</i>	anaerobe	37	WCA
	<i>Collinsella aerofaciens</i>	anaerobe	37	WCA
	<i>Atopobium parvulum</i>	anaerobe	37	PYG
	<i>Eggerthella lenta</i>	anaerobe	37	NB / WCA
	<i>Bifidobacterium longum</i>	anaerobe	37	WCA
	<i>Clostridium ramosum</i>	anaerobe	37	WCA
	<i>Staphylococcus aureus</i>	aerobe	37	TSA/TSB
	<i>Klebsiella pneumoniae</i>	aerobe	28	TSA/TSB
	<i>Escherichia coli</i> LF82	anaerobe	37	WCA
	<i>Shigella flexneri</i>	aerobe	37	NB / TSA
	<i>Oscillibacter valericigenes</i>	anaerobe	30	PYG
	<i>Akkermansia muciniphila</i>	anaerobe	37	Schaedler
	<i>Ruminococcus gnavus</i>	anaerobe	37	WCA
	<i>Bacteroides vulgatus</i>	anaerobe	37	WCA
	<i>Pseudomonas aeruginosa</i>	aerobe	37	TSA/TSB
	<i>Citrobacter freundii</i>	aerobe	37	TSA/TSB
	<i>Enterobacter cloacae</i>	aerobe	28	TSA/TSB
<i>Listeria welshimeri</i>	aerobe	37	TSA/TSB	
<i>Microbacterium flavum</i>	aerobe	28	TSA/TSB	