## **Supplementary Table S1.** Donor Characteristics and Clinical Information

| No. | Sex    | Age, y | Cause of Death           | Time Between Death and Fixation, h |
|-----|--------|--------|--------------------------|------------------------------------|
| 1   | Male   | 52     | Cerebral hernia          | 3                                  |
| 2   | Female | 48     | Intracerebral hemorrhage | 2                                  |
| 3   | Female | 46     | Intracerebral hemorrhage | 1                                  |
| 4   | Male   | 54     | Intracerebral hemorrhage | 2                                  |
| 5   | Male   | 55     | Intracerebral hemorrhage | 2                                  |
| 6   | Male   | 37     | Craniocerebral trauma    | 2                                  |
| 7   | Female | 51     | Craniocerebral trauma    | 2                                  |
| 8   | Male   | 55     | Intracerebral hemorrhage | 1                                  |
| 9   | Female | 64     | Cerebrovascular accident | 3                                  |
| 10  | Female | 60     | Cerebrovascular accident | 3                                  |

## Supplementary Table S2. List of Assays Used in AqH Detection

| Factors  | Assays                     | Company                      | Presented         |
|----------|----------------------------|------------------------------|-------------------|
| ROS      | Chemiluminescence          | Tohoku Electronic Industrial | chemiluminescence |
|          | analyzing system           |                              | counts/10s        |
| 8-OHdG   | Enzyme-linked              | IBL International GmbH       | ng/mL             |
|          | immunosorbent assay        |                              |                   |
| MDA      | Lipid Peroxidation (MDA)   | Bio Vision                   | nmol/L            |
|          | Colorimetric/Fluorometric  |                              |                   |
|          | Assay Kit                  |                              |                   |
| SOD      | SOD Activity Kit           | Nanjing Jiancheng            | U/mL              |
|          |                            | Bioengineering Institute     |                   |
| GSH/GSSG | GSH/GSSG Ratio Detection   | Abcam                        | -                 |
| ratio    | Assay Kit                  |                              |                   |
| SASP     | Multiplex Immunoassay      | Thermo Fisher Scientific     | pg/mL             |
| SA-β-gal | Senescence β-Galactosidase | Cell Signaling Technology    | -                 |
|          | Staining Kit               |                              |                   |

AqH, aqueous humor; 8-OHdG, 8-hydroxydeoxyguanosine; GSH, reduced glutathione; GSSG, oxidized glutathione; MDA, malondialdehyde; ROS, reactive oxygen species; SA-β-gal, senescence-associated beta-galactosidase; SASP, senescence-associated secretory phenotype.

**Supplementary Table S3.** List of Antibodies Used for Immunofluorescence

| Primary   | Type | Company     | Catalog # | Dilution | Secondary   | Company    | Catalog |
|-----------|------|-------------|-----------|----------|-------------|------------|---------|
| Antibody  |      |             |           |          | Antibody    |            | #       |
| 8-OHdG    | M    | Abcam       | ab62623   | 1:200    | Anti-       | Cell       | #4409   |
|           |      |             |           |          | mouse IgG   | Signaling  |         |
|           |      |             |           |          | conjugated  | Technology |         |
|           |      |             |           |          | to Alexa    |            |         |
|           |      |             |           |          | Fluor 555   |            |         |
| p16       | M    | Abcam       | ab108349  | 1:100    | Anti-rabbit | Cell       | #4413   |
|           |      |             |           |          | IgG         | Signaling  |         |
|           |      |             |           |          | conjugated  | Technology |         |
|           |      |             |           |          | to Alexa    |            |         |
|           |      |             |           |          | Fluor 555   |            |         |
| p21       | M    | Novus       | NBP2-     | 1:100    | Anti-       | Cell       | #4409   |
|           |      | Biologicals | 29463     |          | mouse IgG   | Signaling  |         |
|           |      |             |           |          | conjugated  | Technology |         |
|           |      |             |           |          | to Alexa    |            |         |
|           |      |             |           |          | Fluor 555   |            |         |
| p53       | M    | Abcam       | ab237976  | 1:100    | Anti-       | Cell       | #4409   |
|           |      |             |           |          | mouse IgG   | Signaling  |         |
|           |      |             |           |          | conjugated  | Technology |         |
|           |      |             |           |          | to Alexa    |            |         |
|           |      |             |           |          | Fluor 555   |            |         |
| Bcl-2     | M    | Cell        | #3498     | 1:100    | Anti-rabbit | Cell       | #4413   |
|           |      | Signaling   |           |          | IgG         | Signaling  |         |
|           |      | Technology  |           |          | conjugated  | Technology |         |
|           |      |             |           |          | to Alexa    |            |         |
|           |      |             |           |          | Fluor 555   |            |         |
| Bcl-xL    | M    | Novus       | NBP2-     | 1:200    | Anti-       | Cell       | #4409   |
|           |      | Biologicals | 32917     |          | mouse IgG   | Signaling  |         |
|           |      |             |           |          | conjugated  | Technology |         |
|           |      |             |           |          | to Alexa    |            |         |
|           |      |             |           |          | Fluor 555   |            |         |
| Bax       | M    | Invitrogen  | MA5-      | 1:100    | Anti-       | Cell       | #4409   |
|           |      |             | 14003     |          | mouse IgG   | Signaling  |         |
|           |      |             |           |          | conjugated  | Technology |         |
|           |      |             |           |          | to Alexa    |            |         |
|           |      |             |           |          | Fluor 555   |            |         |
| Cleaved   | M    | Cell        | #9664     | 1:500    | Anti-rabbit | Cell       | #4413   |
| caspase-3 |      | Signaling   |           |          | IgG         | Signaling  |         |
|           |      | Technology  |           |          | conjugated  | Technology |         |
|           |      | -           |           |          | to Alexa    | -          |         |
|           |      |             |           |          | Fluor 555   |            |         |
| CREB      | M    | Cell        | #9197     | 1:800    | Anti-rabbit | Cell       | #4413   |

|       |   | Signaling<br>Technology |          |       | IgG<br>conjugated<br>to Alexa | Signaling<br>Technology |       |
|-------|---|-------------------------|----------|-------|-------------------------------|-------------------------|-------|
| SIRT1 | M | Abcam                   | ab110304 | 1:200 | Fluor 555<br>Anti-            | Cell                    | #4409 |
|       |   |                         |          |       | mouse IgG conjugated          | Signaling<br>Technology |       |
|       |   |                         |          |       | to Alexa Fluor 555            | reemieregy              |       |

M, Monoclonal; P, polyclonal.
Bcl-2, B cell lymphoma-2; Bcl-xL, B-cell lymphoma-extra-large; CREB, cAMP response element binding protein; 8-OHdG, 8-hydroxydeoxyguanosine; SIRT1, Sirtuin1.