### Supplementary Methods

## Study design

1

2

3 In this study, volunteers from the Escola Paulista de Medicina/UNIFESP and other institutions in 4 São Paulo were recruited between January and November 2021. The participants were healthy, 5 of both sexes, >18 years old, vaccinated with two homologous doses of CoronaVac (CV, n = 25) 6 or ChAdOx1 (Ch, n = 23), and with a third heterologous dose of BNT162b2 (CV/CV/BNT or 7 Ch/Ch/BNT groups). For CV/CV/BNT, blood samples were collected at five time points: 0-7 days 8 before the first dose (t0), 28 days after the first dose and before the second dose (t1), 14 days 9 after the second dose (t2), 75 days after the second dose (t3), and 14 days after the booster dose 10 with BNT162b2 (t4) (Figure 1A). For Ch/Ch/BNT, blood samples were also collected at five time 11 points: 0-7 days before the first dose (t0), 75 days after the first dose and before the second dose 12 (t1'), 14 days after the second dose (t2'), 75 days after the second dose (t3'), and 14 days after 13 the booster dose with BNT162b2 (t4') (Figure 1A). 14 All participants were tested for anti-S1 IgG at all time points. Neutralizing antibodies and cellular 15 responses were evaluated at t3, t3', t4, and t4'. All individuals enrolled in this study provided written 16 informed consent as part of the protocols approved by the Ethics Committee of the Federal 17 University of São Paulo and by the National Ethics Committee (Comissão Nacional de Ética em 18 Pesquisa) (CONEP, study number CAAE: 32571720.0.0000.5505). The study was conducted in 19 compliance with the principles of the Declaration of Helsinki<sup>1</sup>.

### **Blood collection**

20

21

22

23

24

Venous peripheral blood sample (10 ml) was collected using standard venipuncture and processed within 3 h in the laboratory. Whole blood samples were used immediately for interferongamma (IFN-γ) release assay (IGRA). Plasma samples were obtained after blood centrifugation at 3,000 rpm for 10 min and stored at -80°C for further analysis.

#### **Anti-SARS-CoV-2 antibodies**

Quantitative anti-S1 IgG (Anti-SARS-CoV-2 QuantiVac ELISA IgG) and semiquantitative anti-S1 IgA (Anti-SARS-CoV-2 ELISA IgA), anti-NCP IgG (Anti-SARS-CoV-2 NCP ELISA IgG) levels were measured in plasma using Anti-SARS-CoV-2 ELISA kits (EUROIMMUN, Lübeck, Germany), according to the manufacturer's instructions. Results of quantitative anti-S1 IgG are expressed in binding antibody units/ml (BAU/ml): <25.6 BAU/ml were considered negative, ≥25.6 and <35.2 BAU/ml were indeterminate, and ≥35.2 BAU/ml were positive. The results of semiquantitative anti-S1 IgA and anti-NCP IgG titers are expressed as ratios from the extinction coefficient of the sample to that of the calibrator: <0.8 were considered negative, ≥0.8 and <1.1 were indeterminate, and ≥1.1 were positive.

### **Neutralizing antibodies**

The neutralizing antibody levels were determined in plasma using the SARS-CoV-2 NeutraLISA assay (EUROIMMUN) according to the manufacturer's instructions. In this assay, neutralizing antibodies compete with biotinylated ACE-2 for the SARS-CoV-2 spike S1 domain. Results are expressed as the percentage of inhibition (%IH) of ACE-2 binding to S1: <20% was considered negative, ≥20% and <35% were indeterminate, and ≥35% was positive.

#### IFN-γ assay

Whole blood samples were collected in heparin sample tubes and stimulated with components of the S1 domain (IGRA tube), mitogen (IGRA Stim), or control without activation (IGRA Blank) using the SARS-CoV-2 IGRA Stimulation Tube Set (EUROIMMUN). After 24-h incubation at 37°C ± 1°C, samples were centrifuged at 12,000 g for 10 min, and the plasma was used for the SARS-CoV-2 Interferon Gamma Release Assay (EUROIMMUN). The analyses were performed according to the manufacturer's instructions. Results are expressed in milli-international units per

- 48 milliliter (mIU/mI): <100 mIU/ml were considered negative, ≥100 and <200 mIU/ml were
- 49 indeterminate, and ≥200 mIU/ml were positive.

## Statistical analysis

- 51 Statistical significance tests were performed using the non-parametric Wilcoxon–Mann–Whitney
- 52 test. Correlation coefficients between two variables were quantified using Spearman's rank-order
- correlation. All tests were performed in a two-sided manner using a significance threshold (p) of
- 54 <0.05, with OriginPro 2016 (OriginLab).

55

56

50

# References

- World Medical Association Declaration of Helsinki. *JAMA* 2013;**310**(20):2191. Doi:
- 58 10.1001/jama.2013.281053.