

1 SARS-CoV-2 and other respiratory pathogens are detected in 2 continuous air samples from congregate settings.

3 Mitchell D. Ramuta¹, Christina M. Newman¹, Savannah F. Brakefield¹, Miranda R. Stauss², Roger W.
4 Wiseman^{1,2}, Amanda Kita-Yarbro³, Eli J. O'Connor⁴, Neeti Dahal⁵, Ailam Lim⁵, Keith P. Poulsen⁵, Nasia
5 Safdar⁶, John A. Marx⁶, Molly A. Accola⁶, William M. Rehrauer^{1,6}, Julia A. Zimmer⁷, Manjeet Khubbar⁷,
6 Lucas J. Beversdorf⁷, Emma C. Boehm⁹, David Castañeda⁹, Clayton Rushford⁸, Devon A. Gregory⁸,
7 Joseph D. Yao¹¹, Sanjib Bhattacharyya⁷, Marc C. Johnson⁸, Matthew T. Aliota⁹, Thomas C. Friedrich¹⁰,
8 David H. O'Connor^{1,2*}, Shelby L. O'Connor^{1,2**}

9 ¹ Department of Pathology and Laboratory Medicine, University of Wisconsin-Madison, Madison, WI,
10 USA

11 ² Wisconsin National Primate Research Center, Madison, WI USA

12 ³ Public Health Madison & Dane County, Madison, WI, USA

13 ⁴ EAGLE School, Fitchburg, WI USA

14 ⁵ Wisconsin Veterinary Diagnostic Laboratory, Madison, WI, USA

15 ⁶ University of Wisconsin Hospitals and Clinics, Madison, WI, USA

16 ⁷ City of Milwaukee Health Department Laboratory, Milwaukee, WI

17 ⁸ Department of Molecular Microbiology and Immunology, University of Missouri, School of Medicine,
18 Columbia, MO, USA

19 ⁹ Department of Veterinary and Biomedical Sciences, University of Minnesota, Twin Cities,
20 Minneapolis, MN, USA

21 ¹⁰ Department of Pathobiological Sciences, University of Wisconsin-Madison, Madison, WI, USA

22 ¹¹ Department of Laboratory Medicine and Pathology, Mayo Clinic, Rochester, MN, USA

23 *These authors contributed equally.

24 *Correspondence can be addressed to:

25 Shelby O'Connor

26 slfeinberg@wisc.edu

27 555 Science Drive

28 Madison, WI USA 53711

29 **Abstract**

30 Two years after the emergence of SARS-CoV-2, there is still a need for better ways to assess the risk
31 of transmission in congregate spaces. We deployed active air samplers to monitor the presence of
32 SARS-CoV-2 in real-world settings across communities in the Upper Midwestern states of Wisconsin
33 and Minnesota. Over 29 weeks, we collected 527 air samples from 15 congregate settings and
34 detected 106 SARS-CoV-2 positive samples, demonstrating SARS-CoV-2 can be detected in air
35 collected from daily and weekly sampling intervals. We expanded the utility of air surveillance to
36 test for 40 other respiratory pathogens. Surveillance data revealed differences in timing and location
37 of SARS-CoV-2 and influenza A virus detection in the community. In addition, we obtained SARS-
38 CoV-2 genome sequences from air samples to identify variant lineages. Collectively, this shows air
39 surveillance is a scalable, cost-effective, and high throughput alternative to individual testing for
40 detecting respiratory pathogens in congregate settings.

41 Introduction

42 Viral testing and surveillance have been a challenge throughout the COVID-19 pandemic in the
43 United States. To date, nasal swab-based testing has predominated. Such testing did not reliably and
44 consistently detect SARS-CoV-2 upon its arrival in the United States. First-generation PCR assays
45 were problematic and slow to deploy¹. A massive increase in the need for PCR testing strained
46 supply chains and laboratory capacity, leading to lengthy turnaround times^{2,3}. The development of
47 saliva-based SARS-CoV-2 testing in 2020 reduced the impact of certain supply chain bottlenecks,
48 but laboratory capacity and test availability remained a problem in the United States⁴. Lower-cost,
49 point-of-care antigen tests became available in late 2020 but were not widely used for at-home testing
50 until the arrival of the Omicron Variant of Concern in late 2021. The explosive spread of Omicron was
51 unprecedented; on January 10, 2022, the United States reported the highest single-day COVID-19
52 case number of over 1.3 million new cases and had a seven-day case average three times greater than
53 the previous highest peak recorded in January 2021^{5,6}. Once again, testing laboratories could not scale
54 to meet surging demand and at-home antigen test results, which are rarely reported to public health
55 authorities, were not considered in case counts. As a result, existing testing systems for COVID-19
56 have provided case counts that are, at best, crude estimates for disease burden and transmission risk.
57 Highlighting this discordance, serological data estimates that for every diagnosed case of COVID-19,
58 there are 4.8 undiagnosed SARS-CoV-2 infections⁷.

59 Such swab-based estimates of community infection rates are likely to become less accurate in the
60 future. Antigen testing will cause fewer people to seek out formal diagnostic testing from providers
61 who report test data to public health authorities; indeed, the United States federal government recently
62 began the distribution of 500 million antigen tests without requiring mandatory results reporting⁸. There
63 has been a growing concern that the mental and physical exhaustion caused by COVID-19, often
64 referred to as “pandemic fatigue,” could reduce willingness to seek out testing when symptomatic,
65 particularly if a positive test result is linked to undesirable outcomes such as mandatory isolation⁹.

66 Accurate estimates of SARS-CoV-2 risk are especially important in congregate settings where
67 individuals with varying degrees of risk are in close contact. Highlighting the importance of these
68 settings, the United States invested more than \$12 billion in March 2021 to expand testing in schools,
69 workplaces, long-term care facilities, and underserved congregate settings¹⁰. Evidence of increased
70 SARS-CoV-2 risk from testing programs can be used as an impetus to intensify mitigation measures,

71 such as recommending or requiring facial masking or increased testing. Conversely, such measures
72 can be relaxed when the risk of infection is low. This has led some to advocate for frequent, routine
73 testing of everyone in congregate settings¹¹⁻¹⁴. Many different approaches to high throughput testing
74 have been deployed in support of such comprehensive testing^{4,15-17}, but these are expensive and
75 difficult to scale and maintain. There is also a risk that such resource-intensive testing programs will
76 perpetuate inequalities in the distribution of COVID-19 testing access¹⁸⁻²³.

77 Alternative “environmental” testing strategies that do not rely on individualized testing could provide
78 a more accurate, rapid, and efficient assessment of SARS-CoV-2 infection risk in congregate
79 settings. To date, wastewater-based surveillance for SARS-CoV-2 has received the most attention
80 as an environmental testing strategy²⁴⁻²⁹. Viral RNA is shed in the feces of 30-66% of individuals with
81 COVID-19, regardless of their symptoms³⁰⁻³², allowing SARS-CoV-2 to be detected in wastewater
82 samples. Untreated wastewater collected at municipal wastewater treatment plants includes fecal
83 and liquid waste from households in a sewershed and represents an efficient pooled sample that can
84 provide information on asymptomatic and symptomatic SARS-CoV-2 infections.

85 Wastewater-based surveillance can provide population-wide surveillance data for large geographic
86 regions but mainly relies on fixed sampling locations, limiting its ability to provide spatial resolution
87 within a sewershed. Collecting sewage from individual buildings is possible, and sample collection
88 can be challenging due to differences in the design and complexity of wastewater infrastructure.
89 Furthermore, wastewater surveillance that relies on sample collection from wastewater treatment
90 plants does not capture communities with decentralized systems (e.g., septic tanks) or sites where
91 sewage is pre-treated for decontamination before reaching the wastewater treatment plant (e.g.,
92 hospitals). This is a limitation because the prevalence of COVID-19 and risk of transmission may vary
93 substantially across different congregate settings in a community. There is still a need to develop agile
94 and mobile surveillance technologies to collect hyperlocal data with higher resolution than is possible
95 with wastewater.

96 Air surveillance is an alternative form of environmental sampling for SARS-CoV-2. Passive and active
97 air sampling techniques have been used for air surveillance of viruses, bacteria, and fungi that are
98 released in respiratory droplets and aerosols when infected individuals talk, cough, sneeze, and
99 breathe³³⁻³⁶. Notably, the United States Department of Homeland Security established the BioWatch

100 Program in 2003 to use active air samplers as routine environmental monitoring systems to detect
101 specific biological threats to combat bioterrorism³⁷.

102 Continuous air sampling has key advantages over widespread individual testing and wastewater
103 testing for surveillance in congregate settings. Air surveillance is cost-effective compared to frequent
104 swab-based testing of individuals. Air samples contain a mixture of exhaled components from many
105 individuals and can capture pathogen-containing droplets and aerosols from infectious individuals,
106 enabling virus detection independent of symptoms, test-seeking behavior, and access to swab-based
107 testing. In contrast to wastewater surveillance, active air samplers can be easily moved to different
108 locations, making it possible to collect surveillance data with ultrahigh spatial resolution (e.g., a single
109 room in a building).

110 Several studies have shown the utility of active air samplers to detect aerosols containing SARS-
111 CoV-2³⁸⁻⁴³ in controlled settings and locations with known SARS-CoV-2 cases. Horve et al. (2021)
112 demonstrated consistent detection of heat-inactivated SARS-CoV-2 virus at an aerosol concentration
113 of 0.089 genome copies per liter of air (gc/L) when air samples were collected in a room-scale
114 experiment during an eight-hour interval³⁸. Another study compared the effectiveness of surface and
115 bioaerosol sampling methods to detect SARS-CoV-2 and showed active air samples detected SARS-
116 CoV-2 in 53% of the samples when run for 1-2 hours in hospital rooms of COVID-19 patients, while
117 passive air sampling and surface swabs detected SARS-CoV-2 in only 12% and 14% of samples,
118 respectively³⁹. Parhizkar et al. (2021) used active air samplers to assess the relationship between
119 COVID-19 patient viral load and environmental viral load in a controlled chamber over three days.
120 Increases in patient viral load were associated with lower cycle threshold (Ct) values detected in near
121 (1.2 meters) and far (3.5 meters) air samplers⁴⁰. Lastly, a study demonstrated the utility of using active
122 air samplers to track the presence and concentration of virus in air longitudinally during COVID-19
123 isolation periods in student dormitories. The study observed a significant increase in Ct values for
124 COVID-19 positive students after their first test, as well as in environmental samples as individuals
125 recovered indicating a reduction in virus presence⁴¹. These studies provide proof of concept on the
126 feasibility of using active air samplers to detect SARS-CoV-2. However, each study was performed in a
127 controlled environment occupied by COVID-19 positive individuals.

128 Here, we evaluate whether active air samplers can be used for prospective air surveillance of SARS-
129 CoV-2 in real-world congregate settings, where pathogen-containing aerosols are likely present at a

130 much lower concentration, and the presence of positive individuals is unknown. This study addresses
131 a key knowledge gap for how active air samplers perform as routine pathogen monitoring systems
132 in real-world settings. We demonstrate that it is feasible to use active air samplers for air respiratory
133 pathogen detection and sequencing across different types of congregate settings.

134

135 **Results**

136 Study design

137 From July 19, 2021, to February 9, 2022, continuous air samplers were deployed in several public
138 locations to survey SARS-CoV-2 in the environment of real-world settings. We used Thermo Fisher
139 Scientific AerosolSense Samplers for daily and weekly air surveillance at places considered to be
140 high-risk for close-contact, indoor SARS-CoV-2 transmission. Air cartridges were collected and
141 tested for SARS-CoV-2 RNA by quantitative reverse transcription PCR (RT-qPCR) and transcription-
142 mediated amplification (TMA) (Figure 1). Several different RT-qPCR and TMA assays were used to
143 test samples for viral RNA (vRNA) throughout the study, depending on the location of the test site and
144 the availability of testing at the time. Further details on the cut-off values used for calling air samples
145 positive, inconclusive, or negative are described in the methods section. We developed a user-friendly
146 workflow to collect air cartridge metadata, upload test data, and report results to surveillance sites
147 within 24-48 hours of receiving air cartridges for testing (Figure 1), enabling non-technical staff (e.g.,
148 custodial staff, students) to exchange and catalog cartridges easily and accurately.

149

150 SARS-CoV-2 detection in community settings

151 To demonstrate the utility of air surveillance in real-world settings, we chose a diversity of community
152 locations for placement: a campus coffee shop, hospital, office, campus athletic training facility,
153 brewery taproom, cafeteria, bar, two preschools, four K-12 schools, and two shelters located
154 throughout Wisconsin and Minnesota. We collected a total of 527 air cartridges from the 15 locations
155 to test for the presence of SARS-CoV-2 genomic material (Table 1). Four hundred sixty-six (88.4%)
156 air samples were collected from testing sites in Dane County, Wisconsin, 26 (4.9%) from Minneapolis,
157 Minnesota, 19 (3.6%) from Rochester, Minnesota, and 16 (3.1%) from Milwaukee, Wisconsin. During
158 the 29 weeks, Dane County experienced a moderate-to-high transmission rate of COVID-19 cases

159 despite having a high county-wide vaccination rate (>65% adults with two doses) (Supplementary
160 Figure 1). Increases in COVID-19 cases were observed following the emergence of Delta and Omicron
161 Variants of Concern in the community. An emergency mask mandate was instituted by Public
162 Health Madison Dane County on August 19th, 2021 and extended throughout the entire study. The
163 order required every individual aged two and older to wear a face-covering in most public enclosed
164 spaces, including K-12 schools. Exceptions were made in spaces where all people were known to be
165 vaccinated.

166 Throughout the study, we detected a total of 106 SARS-CoV-2 positive and 52 inconclusive air
167 cartridges (an inconclusive result is defined when at least one of the PCR targets is positive while at
168 least one of the PCR targets is negative). We were able to intermittently identify SARS-CoV-2 positive
169 air cartridges at 14 of the 15 surveillance sites, even when intensive risk mitigation strategies were
170 implemented by public health (Table 1). Together, this provides evidence that SARS-CoV-2 can be
171 detected in continuous air samples collected from a variety of prospective, real-world settings.

172 By testing air samplers in these real-world settings, it was impossible to know the SARS-CoV-2 status
173 of every person who spent time in the proximity of the samplers, but we were able to retrospectively
174 correlate air surveillance data with reported test results during a prolonged COVID-19 outbreak at one
175 of the testing sites (Figure 2). The first confirmed case was reported by an individual who quarantined
176 at home after experiencing symptoms on day 10, had a diagnostic test specimen collected two
177 days later, and received test results on day 14. Individuals who congregated in the same room were
178 quarantined after the case was reported and no one from the room tested positive during that time.
179 Three colleagues, considered to be close contacts with the first COVID-19 positive person but worked
180 in different rooms, tested positive with antigen tests while in the building on days 14, 15, and 16. Two
181 of the individuals were symptomatic and one was asymptomatic at the time of testing. Two more
182 rooms were quarantined and seven people tested positive while at home. A fifth individual tested
183 positive on-site on day 23 after being in the building for the entire prior week; close contacts in their
184 room were quarantined and eight individuals tested positive while in quarantine. Air samples were
185 positive before the first confirmed case of COVID-19 and throughout the outbreak that resulted in a
186 total of 20 confirmed cases; SARS-CoV-2 genetic material was detected in an air sample for the first
187 time between days 5 and 7, preceding the first confirmed case by 7 days. Air samples collected after
188 day 23 were either inconclusive or negative for SARS-CoV-2. No reported cases were observed in the

189 building during this time. No air sample was collected from day 28-30 because the air cartridge was
190 inserted improperly, leading to a machine error. This suggests that congregate risk during this outbreak
191 could have been estimated with air sampling data alone had individual testing not been available or
192 widespread, and that air sampling provided an early indication of SARS-CoV-2 transmission risk. It
193 should be noted that, at the time of this outbreak, we did not have sufficient data on air sampling
194 accuracy to make recommendations to the affected congregate setting to intensify risk mitigation.
195

196 Extending the duration of continuous air sampling

197 Daily air surveillance programs are resource-intensive and expensive. To reduce the cost and
198 complexity of surveillance programs, we tested the feasibility of extending the sampling interval.
199 Over the course of five weeks, two adjacent AerosolSense instruments were deployed to either
200 run continuously (~168 hours) or daily (~24 hours) over several days. Air samples were gathered
201 throughout the week, nucleic acids were isolated simultaneously, and tested with two RT-qPCR Center
202 for Disease Control and Prevention (CDC) assays targeting SARS-CoV-2 N1 and N2. During the first
203 week of sampling, both air samplers detected SARS-CoV-2 genetic material, showing continuous air
204 sampling for a week captured similar SARS-CoV-2 genetic material, determined by Ct values, when
205 compared to more resource-intensive daily sampling (Figure 3). Weekly and daily air sample results
206 were concordant for the next two weeks; air samples were either negative or inconclusive. Continuous
207 air samples were negative for the last two weeks, while daily sampling identified one inconclusive
208 sample during each of these weeks. However, each of the inconclusive samples had high N2 Ct
209 values of 39.1 and 39.4, respectively. This suggests that congregate settings can use either daily or
210 weekly sampling to achieve a balance between cost and turnaround time while maintaining detection
211 sensitivity.

212

213 Expanding detection to additional respiratory pathogens

214 To explore whether pathogens other than SARS-CoV-2 are detected in the same collected air sample,
215 we tested them for the presence of nucleic acids of 40 other viral, bacterial, and fungal respiratory
216 pathogens using the TrueMark Respiratory Panel 2.0 Array Card on a QuantStudio 7 Pro Instrument⁴⁴.
217 We assessed the limit-of-detection (LOD) for each TrueMark assay using contrived air samples to

218 determine pathogen-specific thresholds for calling samples positive or negative. Contrived air samples
219 were created and tested in quadruplicate by spiking the TrueMark Respiratory Panel 2.0 Amplification
220 Control (Thermo Fisher Scientific) plasmid into pooled air samples at concentrations of 0, 1.25, 12.5,
221 50, and 250 copies per reaction. Air samples were collected for 48 hours from an empty office to
222 minimize any background pathogens present in the samples (see methods section for more details).
223 Pathogen-specific thresholds were determined by averaging the Cycle relative threshold (C_{rt}) values
224 of positive replicates at the lowest dilution concentration with at least 75% positive replicates. Cycle
225 relative threshold values listed in Supplementary Table 1 were used as cut-off values for positivity for
226 the detection of each pathogen.

227 From October 25, 2021, to February 9, 2022, air samples were collected weekly from eight sites in
228 Dane County, Wisconsin: a campus coffee shop, a preschool, an office space, a campus athletic
229 training facility, and four K-12 schools. Semi-quantitative RT-PCR was performed on 104 air samples
230 using this TrueMark Respiratory 2.0 Panel. During the 15 weeks of air surveillance, we detected
231 16 different respiratory pathogens across the eight sites (Supplementary Table 2). Commensal or
232 transiently commensal respiratory tract microbes were frequently detected in air samples at each
233 site, including *Klebsiella pneumoniae*, *Staphylococcus aureus*, and *Moraxella catarrhalis* (Figure
234 4). The panel also detected respiratory pathogens associated with illness in school-aged children,
235 including adenovirus, human coronavirus OC43, Epstein-Barr virus, cytomegalovirus, influenza A
236 virus, parainfluenza virus 3, respiratory syncytial virus A, and SARS-CoV-2. Certain pathogens, such
237 as human bocavirus, were consistently detected only in settings where there were young children,
238 consistent with its widespread distribution in this population, highlighting that different types of
239 congregate settings have distinctive air surveillance pathogen signatures (Figure 4B)⁴⁵.

240 The pattern of influenza A virus (IAV) nucleic acid detection was especially striking. The Wisconsin
241 Department of Health Services issued a health alert on November 30, 2021, noting increased IAV
242 activity among college and university students⁴⁶. As shown in Figure 5, we detected IAV nucleic acids
243 in the air collected from two AerosolSense samplers deployed on the University of Wisconsin-Madison
244 campus beginning the week of November 10, 2021. Air collected from both of these samplers
245 contained IAV nucleic acids from mid-November 2021 to January 2022. One sampler, located in
246 a campus coffee shop, was negative for IAV beginning the week of December 22, 2021, through
247 the week ending January 12, 2022 (Figure 5A); this coincided with the end of the UW-Madison fall

248 academic semester and holiday break. During this time, the coffee shop was closed for customers
249 from December 17, 2021, to January 18, 2022, but the building was still open for repairs and
250 cleaning. We continued to detect IAV nucleic acids collected by the sampler in the training facility, as
251 student-athletes and staff remained on campus for training and competition during the holiday break.
252 Strikingly, IAV nucleic acids were only sporadically detected in air samplers located on non-campus
253 community sites in Dane County, Wisconsin. In contrast, SARS-CoV-2 nucleic acids were frequently
254 detected at testing sites across the community (Figure 5B). Dane County experienced a high
255 transmission rate of COVID-19 during this time and the detection of SARS-CoV-2 positive air samples
256 increased at testing sites following the emergence of Omicron. Overall the differential detection of IAV
257 and SARS-CoV-2 nucleic acids is consistent with the localized known IAV campus outbreak and wide-
258 spread SARS-CoV-2 transmission.

259

260 Sequencing of SARS-CoV-2 RNA from collected air samples

261 Throughout the pandemic, deep sequencing of SARS-CoV-2 RNA extracted from clinical samples and
262 wastewater has played a crucial role in monitoring viral evolution and tracking variants of concern.
263 We used two sequencing strategies to obtain partial and near-full genome SARS-CoV-2 sequences
264 from 11 air samples to demonstrate the feasibility of genotyping SARS-CoV-2 from collected air.
265 Sequencing efforts were focused on air samples with low Ct values from congregate settings where
266 individuals often removed their masks to eat and drink (e.g. taproom and shelters).

267 Targeted sequencing was performed on nine samples collected from two AerosolSense samplers in
268 a brewery taproom in Minnesota between November 22, 2021, and January 25, 2022, using primers
269 targeting the SARS-CoV-2 spike gene receptor-binding domain (RBD)²⁶. Consensus sequences from
270 four samples collected between November 22 and December 13, 2021, all contained the characteristic
271 S:L452R and S:T478K variants associated with the Delta lineage that predominated at this time
272 (Table 2). Interestingly, one cartridge collected genetic material that had two rare consensus variants,
273 S:F456L and synonymous S:F562F, suggesting that sequencing can detect unexpected variants
274 in air samples and could be used for detecting newly emerging variants of concern in congregate
275 settings. Six air samples collected between December 30, 2021, and January 25, 2022, all contained
276 characteristic S:K417N, S:N440K, S:G446S, S:S477N, S:T478K, S:E484A, S:Q493R, S:G496S,
277 S:Q498R, S:N501Y, S:Y505H, S:T547K variants associated with the Omicron BA.1 lineage coinciding

278 with the emergence of Omicron in the region (Table 2). In early January, we detected both Delta and
279 Omicron sequences in one of the air samples. These data support that virus genetic material collected
280 by air samples parallels the longitudinal detection of the same lineages transmitted in the community.
281 The remaining two additional samples with Ct values below 32 collected from two shelters in
282 Milwaukee, Wisconsin, from December 21, 2021, to January 7, 2022, were examined using the ARTIC
283 protocol followed by Illumina sequencing. Sequence coverage was incomplete for both samples, with
284 28% and 4% of the sequences having low coverage or missing data. However, there was enough
285 information to assign both samples to the Omicron BA.1 lineage (Table 2).
286

287 **Discussion**

288 An increasing number of public health organizations have employed environmental surveillance
289 methods, in conjunction with individual clinical testing, to provide data on SARS-CoV-2 prevalence,
290 monitor viral evolution, and track viral variants in communities^{24,29,47,48}. Environmental surveillance tools
291 could help public health make data-driven decisions for implementing COVID-19 interventions and
292 allocating resources within a community. This study demonstrates the feasibility of using active air
293 samplers for environmental pathogen surveillance in real-world settings, even while intensive risk-
294 mitigation strategies are implemented. Air surveillance may improve our ability to identify pathogen-
295 containing aerosols present in congregate settings to assess the risk of transmission with high spatial
296 resolution, providing a more complete picture to public health officials.

297 An especially attractive element to air surveillance is its relative cost-effectiveness compared to
298 intensive, individual testing in congregate settings. For example, Wisconsin received more than
299 \$175 million in United States federal funding to support SARS-CoV-2 testing of individual teachers,
300 staff, and students in 2021⁴⁹. There are 5,987 public and private elementary and secondary schools
301 in Wisconsin, plus 2,716 preschool and child development centers^{50,51}. These are geographically
302 distributed throughout the state, such that air sampling results in schools could be generally
303 representative of the communities in which they are located. At approximately \$5,000 per air sampler
304 (<https://www.thermofisher.com/order/catalog/product/AEROSOLSENSE>), each school in the state
305 could be equipped with an air sampler for approximately \$43.5 million. Assuming daily testing for
306 180 school days per academic year at a cost of \$115 per sample (\$40 per cartridge and \$75 per

307 test), a testing program would be cost-similar to the investment in individual testing (\$180 million)⁵².
308 A more strategic design that tests daily during the high-risk period for respiratory pathogens (ie.,
309 mid-November through February) and weekly during the remainder of the school year would be
310 considerably less expensive. Moreover, the possibility of substituting low-cost nylon flocking material
311 for the AerosolSense Capture Material, as described by others, could enable cartridge recycling,
312 substantially reducing both cost and plastic waste³⁸. While there would undoubtedly be other
313 unforeseen costs of deploying air sampling at scale (e.g., creation of high volume data management
314 solutions), the ability to task-shift air sampling cartridge management to staff with no scientific training
315 makes air sampling more accessible and scalable than individual testing with specialized personnel.
316 Furthermore, air sampling programs could help lessen the burden on individuals participating in testing
317 programs compared to current swab-based surveillance testing.

318 Networks of air samplers deployed as described here could also play a role in improving public
319 health resilience to new and emerging respiratory diseases. Had a nationwide network existed in
320 the United States prior to the arrival of SARS-CoV-2, the spread of the virus through space and time
321 could have been more accurately evaluated. Moreover, adding continuous air sampling in the future
322 at ports of entry and, potentially, aboard international aircraft could intercept vessels and passengers
323 harboring worrisome respiratory pathogens. Improvements in technologies that enable real-time,
324 highly multiplexed pathogen detection and genotyping could be leveraged with air sampling to
325 improve quarantine effectiveness. Consider the arrival of the Omicron variant: it was initially reported
326 to the World Health Organization by South African authorities on November 24, 2021⁵³. The first
327 confirmed Omicron case, however, occurred weeks earlier in a sample collected on November
328 9, 2021. Additionally, continuous genomic surveillance of air by targeted spike or whole-genome
329 sequencing from international travelers arriving in the United States could have shortened the window
330 of first detection of Omicron in this country. Establishing a network of air surveillance programs now
331 could provide an early warning for the arrival of future SARS-CoV-2 variants, as well as future novel
332 respiratory viruses of concern.

333 We used RT-qPCR analysis for detecting low levels of pathogen nucleic acids captured in air samples.
334 Quantitative RT-PCR assays of nucleic acids in a test tube provide semi-quantitative information on
335 the viral load present in ambient air, but these viral load data should be interpreted with care. The
336 exact environmental viral load and duration of exposure are unknown because of the large sampled air

337 volumes collected during daily and weekly sampling intervals in this study. Furthermore, the amount
338 of genetic material collected by active air samplers depends on many factors, including the location of
339 the sampler, ventilation capacity, the amount of virus shed by each infected individual, the number of
340 infected individuals in an area, and the dimensions of the indoor environment. We are now optimizing
341 the placement of air samplers by evaluating carbon dioxide levels, air exchange, and movement
342 data with low-cost sensor arrays to define locations most likely to sample air from a large number of
343 individuals in a congregate setting.

344 Additionally, RT-qPCR analysis cannot distinguish whether an air cartridge is collecting infectious virus
345 or genetic material that does not pose a risk for infection. We did not attempt to culture SARS-CoV-2
346 or any other pathogens from air cartridges in this study to determine if infectious virus can be isolated
347 from the nucleic acid positive samples. Several studies have attempted to culture SARS-CoV-2 and
348 other viruses from air samples with mixed results^{42,43,54,55}. Furthermore, previous studies have shown
349 that the air sampling method and capture media can affect viral integrity⁵⁶. AerosolSense samplers use
350 an accelerating slit impactor to collect aerosol particles on dry filter capture substrate. Air sampling
351 methods that rely on impactors and filters are not optimal for maintaining virus viability because
352 of damage caused by impaction forces and dehydration during the collection process. Live virus
353 recovery from continuous air samples would be valuable, as it might potentiate culture and expansion
354 of unknown pathogens with greater sensitivity.

355 After two years of COVID restrictions, there is pushback against public health measures to counteract
356 virus transmission^{9,57}. One component of this resistance is that guidelines issued at the national,
357 state, or even municipal level do not necessarily reflect hyperlocal risk within specific congregate
358 settings: an individual school, sports arena, bar, etc. Air sampling provides a measure of risk with this
359 level of granularity. However, care must be taken when interpreting and sharing air sampling data. In
360 some settings, stakeholders may choose to be liberal in disclosing air sampling results, sharing this
361 information with employees, customers, visitors, and others so they can individually modulate their
362 risk mitigation. In other settings, public health and testing laboratories may work directly with the
363 leadership of congregate settings to couple air sampling data with action. For example, one of the
364 county public health departments involved in this work offers enhanced testing to sites where high
365 levels of SARS-CoV-2 is detected in the air, while a second public health department created a flow

366 chart describing how schools might respond to positive air sampling data if there are no known cases
367 of SARS-CoV-2 in a given school (Supplementary Figure 2).

368 Additionally, the brewery taproom implemented a vaccine mandate on December 10, 2021, after
369 observing an increase in SARS-CoV-2 detection in air samples and COVID-19 cases in the community.
370 The vaccine mandate required customers to show proof that they received all recommended doses in
371 their primary series of COVID-19 vaccines or had a negative COVID-19 test within the last 72 hours for
372 indoor dining. No air samples collected between December 13, 2021, and December 30th, 2021, were
373 positive for SARS-CoV-2. However, SARS-CoV-2 was detected in 72% of the air samples collected
374 following the emergence of Omicron. These data reiterate that vaccinated individuals infected with
375 SARS-CoV-2 can shed virus⁵⁸.

376 We did not have sequencing data available from the brewery taproom in real-time. In retrospect,
377 sequencing and RT-qPCR data could have been used to help make data-driven decisions to adapt
378 the risk mitigation strategy. Adjustments to the COVID-19 policy could have included increasing the
379 ventilation or expanding the vaccine mandate to require booster doses that have been shown to
380 improve protection against Omicron^{41,59}. However, even with these data available, congregate settings
381 may be hesitant to increase risk mitigation strategies past the most stringent guidelines set out by
382 the CDC and local public health. Furthermore, some settings may have no appetite for COVID-19
383 risk mitigation regardless of air surveillance results. Environmental surveillance in these settings may
384 nonetheless be valuable to public health alone, allowing them to anticipate and respond quickly
385 to surges in respiratory disease⁶⁰. In fact, in settings where diagnostic testing for SARS-CoV-2 is
386 limited by pandemic fatigue and apathy towards risk minimization measures, air sampling could be
387 exceptionally useful in providing baseline data on respiratory virus levels that would otherwise be
388 impossible to obtain.

389 Taken together, these results show that continuous air surveillance with active air samplers can
390 unambiguously detect respiratory pathogens, including SARS-CoV-2, in congregate settings. Similar
391 to the National Wastewater Surveillance System recently established by the US Centers for Disease
392 Control and Prevention, expansion of air surveillance efforts could provide additional safeguards for
393 congregate settings and improve resilience to future respiratory virus threats²⁴.

394 **Methods**

395 Collection of air samples

396 AerosolSense instruments (Thermo Fisher Scientific, cat. 2900AA) were deployed in various indoor
397 community settings for air pathogen surveillance. Samplers were placed in high-traffic areas on flat
398 surfaces 2-5 feet off the ground and calibrated to sample 200 liters of air per minute. AerosolSense
399 cartridges (Thermo Fisher Scientific, cat. 12148001) were installed and removed from the air sampler
400 according to the manufacturer brochure and transferred to the lab for testing⁶¹. We developed a
401 workflow to simplify data collection, management, and reporting (Figure 1). The workflow relies on
402 the iOS and Android Askidd mobile app to easily collect air cartridge metadata and upload it to a
403 centralized LabKey database. Air sampler users simply open the Askidd app, take a picture of the
404 air cartridge barcode when installed and removed from the machine. The Askidd app collects GPS
405 coordinates of the air sampler, timestamp, AerosolSense instrument ID, and air cartridge barcode
406 to send to LabKey. When air sample testing was completed in the lab the results were uploaded to
407 the Labkey database and displayed in the Askidd mobile app. This workflow tracks data for every
408 cartridge and limits user errors that could occur during manual input.

409

410 Detection of SARS-CoV-2

411 University of Wisconsin-Madison

412 Hologic Aptima SARS-CoV-2 Assay

413 AerosolSense cartridges collected at a hospital in Dane County, Wisconsin from August 20th to
414 October 25th, 2021 were tested for SARS-CoV-2 viral RNA using the Aptima SARS-CoV-2 Assay
415 (Hologic) on the Panther System (Hologic). The Aptima SARS-CoV-2 Assay was authorized for
416 emergency use authorization (EUA) by the United States Food and Drug Administration (FDA) for
417 the qualitative detection of vRNA⁶². Air cartridges were collected from AerosolSense Samplers as
418 recommended by the manufacturer. One substrate was removed from the cartridge using sterile
419 forceps, transferred to a tube containing 750 μ L of universal transport medium (Copan), and incubated
420 at room temperature for 5-10 minutes. Following the incubation, 500 μ L of the eluate was transferred
421 to a Panther Fusion Specimen Lysis Tube (Hologic) containing 710 μ L of specimen transport medium.
422 The tube was gently mixed by inverting it several times before loading it into the Panther System to

423 automatically run the Aptima SARS-CoV-2 Assay as described by the manufacturer. Aptima SARS-
424 CoV-2 positive and negative controls were run with each set of air samples. An estimated cut-off value
425 of >650 RLU was used to consider samples as SARS-CoV-2 positive.

426

427 *CDC SARS-CoV-2 RT-qPCR Assay*

428 RNA extraction and real-time reverse-transcriptase polymerase chain reaction (RT-qPCR) testing
429 of the air filter specimens occurred at the University of Wisconsin-Madison WVDL-WSLH COVID
430 Laboratory (WWCL, Madison, Wisconsin). Each air cartridge was submerged in 500 μ L of 1X PBS
431 for at least one hour. For all air cartridges, the tubes were vortexed vigorously and 190 μ L of the PBS
432 was used for RNA extraction using the MagMAX Viral/Pathogen II (MVP II) Nucleic Acid Isolation Kit
433 (Thermo Fisher Scientific) on a 96-well KingFisher Flex extraction platform and eluted in a volume
434 of 50 μ L according to manufacturer's instructions. A multiplex one-step RT-qPCR assay targeting
435 the 2019-nCoV N gene sequences (N1 and N2) and the human RnaseP (RP) gene was used for the
436 SAR-CoV-2 viral detection. The RT-qPCR primers and probes sequences were based on the CDC
437 assay with alternative fluorophores on the 5' end of each probe (along with 3' black hole quencher)
438 for multiplexing⁶³. The N1 probe was labeled with ABY dye, the N2 probe with the FAM dye, and the
439 RP probe with the VIC dye. The 16 μ L reaction mix consists of 1x TaqPath 1-Step Multiplex Master
440 Mix (Thermo Fisher Scientific), 250 nM of each forward and reverse primers for the N1 and N2 targets,
441 100 nM of each forward and reverse primer for the RP gene target, 62.5 nM for each of the N1 and N2
442 probes, 50 nM of the RP probe and 5 μ L of sample RNA or controls. The RT-qPCR amplification was
443 performed with one cycle at 53°C for 10 mins and 95°C for 2 mins, followed by 40 cycles of 95°C for 3
444 secs and 60°C for 30 secs on a QuantStudio 7 Pro Real-Time PCR System (Thermo Fisher Scientific,
445 Inc.). The data was analyzed in the Design and Analysis v2.4 software (Thermo Fisher Scientific, Inc.)
446 using the auto baseline and threshold settings at 0.15 for N1 and N2 and 0.1 for the RP. Samples with
447 amplification (Ct <40) in both the N1 and N2 targets were determined as positive for SARS-CoV-2.
448 In contrast, samples with amplification in only 1 of the targets were determined as inconclusive, and
449 samples without amplification in both N1 and N2 targets were deemed negative for SARS-CoV-2.
450 Each run included a negative extraction control using a pool of previously identified SARS-CoV-2
451 negative samples, positive extraction control, negative template control, and positive amplification
452 control plasmid. The RP gene was utilized for analysis of human nasal swab samples performed at the

453 same laboratory. It was not factored in for the result criteria for the air filter samples due to the low and
454 inconsistent level of human cellular material trapped by the air filters.

455

456 *TrueMark Respiratory Panel*

457 AerosolSense cartridges collected from community testing sites from October 25, 2021 to February
458 9, 2022 were tested for the presence of 41 different respiratory tract viral, bacterial, and fungal nucleic
459 acids using the TrueMark Respiratory Panel 2.0 TaqMan Array Card (TAC) (Thermo Fisher Scientific,
460 Inc.). Substrates were extracted from the AerosolSense cartridge using sterile forceps, submerged
461 into tubes containing 500 μ L of PBS, vortexed for 5 seconds, and stored at 4C for 10-30 minutes.
462 Samples were removed from 4C and sterile forceps were used to disrupt the substrate by pressing
463 it against the bottom of the tube several times to ensure bound particles were eluted into the PBS.
464 According to the manufacturer's recommendations, the substrate was removed from the tube and
465 nucleic acids were isolated from the eluate using the Maxwell Viral Total Nucleic Acid Purification Kit
466 (Promega) with the Maxwell 16 instrument (Promega). Briefly, 300 μ L of the eluate was transferred
467 to a tube containing 300 μ L of lysis buffer and 30 μ L of Proteinase K. A nuclease-free water control
468 was processed with each Maxwell run and used in the TrueMark protocol as a no-template control.
469 Tubes were vortexed for 5 seconds and incubated on a heat block at 56°C for 10 minutes. Following
470 incubation, samples were centrifuged for 1 minute to pellet any debris. Then 630 μ L of each reaction
471 mix was transferred into a Maxwell 16 cartridge, loaded into Maxwell 16 instrument, and processed
472 with the Viral Total Nucleic Acid program. Nucleic acids were eluted in 50 μ L of RNase-free water. To
473 perform the preamplification protocol, 5 μ L of isolated nucleic acids were transferred into a PCR strip
474 tube containing 2.5 μ L of TaqPath 1-Step RT-qPCR Master Mix, CG (Thermo Fisher Scientific, Inc.),
475 and 2.5 μ L of TrueMark Respiratory Panel 2.0 PreAmp Primers (Thermo Fisher Scientific, Inc.). Pre
476 Amplification was performed on a thermocycler with the following cycling conditions: UNG incubation
477 step at 25°C for 2 minutes, reverse transcription at 50°C for 30 minutes, UNG inactivation at 95°C
478 for 2 minutes, 14 cycles at 95°C for 15 seconds (denaturation), 60°C for 2 minutes (annealing and
479 extension), followed by inactivation at 99.9°C for 10 minutes, and 4°C until samples were ready for
480 use. Preamplified products were diluted 1:20 in nuclease-free water, and the TrueMark Respiratory
481 Panel 2.0 Amplification Control (Thermo Fisher Scientific, Inc.) was diluted 1:2 to include with every
482 set of samples. TrueMark reaction mix was prepared by combining 20 μ L of each diluted preamplified

483 product with 50 μ L of TaqMan Fast Advanced Master Mix (Thermo Fisher Scientific, Inc.) and 30 μ L of
484 nuclease-free water. TAC were equilibrated to room temperature, and 100 μ L of each reaction mix was
485 loaded into its respective TAC port. TAC were centrifuged twice at 1,200 rpm for 1 minute each spin.
486 TAC were sealed with a TAC Sealer, loaded into the QuantStudio 7 Pro (QS7), and run with the settings
487 recommended by the manufacturer. Data were exported from the QS7 into the Thermo Fisher Design
488 and Analysis Software 2.6.0. Data were analyzed according to the manufacturer's recommendations
489 using the relative quantification module with the relative threshold algorithm (Ct). Results were
490 exported from the quality check module. Analysis was performed using a custom R script (v. 3.6.0)
491 in RStudio (v. 1.3.959) to filter amplified results using the following cut-off values: amplification score
492 >1.2 and Ct confidence >0.7 . Samples were further filtered on reaction-specific Ct cut-off values
493 determined in a limit-of-detection experiment (Supplementary Data 1). The TrueMark Respiratory Panel
494 2.0 includes technical control assays for human RNase P (RPPH1) and human 18S ribosomal RNA.
495

496 *TrueMark Respiratory Panel Limit-of-Detection Estimation using Contrived Air Samples*

497 Contrived air samples were prepared using the TrueMark Respiratory Panel 2.0 Amplification
498 Control (Thermo Fisher Scientific, Inc.) to estimate the limit of detection (LOD) of the TrueMark
499 Respiratory Panel 2.0 TaqMan Array Card. Briefly, two air cartridges were collected in an office for
500 48 hours each. The air cartridge substrates were processed and total nucleic acids were isolated as
501 described above. Eluates from the four substrates were pooled together and aliquoted into five tubes.
502 TrueMark amplification control plasmid, initially diluted in nuclease-free water, was added to four
503 tubes at dilutions of 50 copies/ μ L, 10 copies/ μ L, 2.5 copies/ μ L, and 0.25 copies/ μ L. Final template
504 concentrations for the preamplification reaction were 250, 50, 12.5, and 1.25 copies per reaction,
505 respectively. No amplification control plasmid was added to the fifth tube that was used to determine
506 the targets present in the background of contrived air samples collected from the empty office.
507 An unused air cartridge was processed with the air samples as a negative template control. Four
508 replicates of each contrived sample and control were processed through the reverse transcription, pre-
509 amplification, dilution, and PCR protocols as described above. Data were analyzed in Thermo Fisher
510 Design and Analysis Software 2.6.0 according to the manufacturer's recommendations. Replicates
511 were called positive using the following cut-off values: amplification score >1.2 and Ct confidence
512 >0.7 . Cycle relative threshold cut-off values were determined by averaging the Ct values of positive

513 replicates at the lowest dilution concentration with at least 75% positive replicates. Any reaction
514 targets that were detected in the contrived air sample or not detected at the highest amplification
515 control dilution were excluded from the analysis and Crt cut-off values defaulted to the manufacturer's
516 recommendation of Crt >30.

517

518 University of Minnesota

519 Nucleic Acid Extraction and RT-qPCR

520 To elute the sample from the AerosolSense cartridges, both substrates were placed into 1mL of
521 PBS, making sure the substrates were fully saturated with PBS. A pipette was used to push down
522 on the substrates to extract as much eluate out of them as possible. Eluate was then transferred to
523 a new tube. Samples were extracted using the *Quick-RNA Viral Kit* (Zymo Research). The extraction
524 method followed manufacturer-recommended protocols with the notable exceptions of using 100 μ L
525 of starting material and eluting with 65 μ L of appropriate elution material as indicated by manufacturer
526 protocols. RT-qPCR reactions were set up in a 96-well Barcoded plate (Thermo Scientific) for either
527 the N1 or N2 primers and probes with CDC-recommended sequences⁶⁴. Then 5 μ L extracted RNA
528 was added to 15 μ L qPCR master mix comprised of the following components: 8.5 μ L nuclease-free
529 water, 5 μ L TaqMan™ Fast Virus 1-Step Master Mix (Thermo), and 1.5 μ L primer/probe sets for either
530 N1 or N2 (IDT, Cat# 10006713). SARS-CoV-2 RNA Control was obtained from Twist Biosciences
531 (Genbank Ref. No. MN908947.3) and used as a positive control in each run. Reactions were cycled in
532 a QuantStudio QS3 (ThermoFisher) for one cycle of 50°C for 5 minutes, followed by one cycle of 95°C
533 for 20 seconds, followed by 50 cycles of 95°C for 3 seconds and 55°C for 30 seconds. A minimum
534 of two no-template controls (NTCs) were included on all runs. Baselines were allowed to calculate
535 automatically, and a ΔR_n threshold of 0.5 was selected and set uniformly for all runs. Ct values were
536 exported and analyzed in Microsoft Excel. Amplification curves were manually reviewed. Samples with
537 Ct<40 in both N1 and N2 reactions were determined as positive for SARS-CoV-2. In contrast, samples
538 with Ct<40 in only N1 or N2 targets were determined as inconclusive results, and samples without
539 amplification in both N1 and N2 targets were deemed negative for SARS-CoV-2.

540

541 *SARS-CoV-2 Spike Receptor Binding Domain Sequencing*

542 Targeted sequencing of SARS-CoV-2 spike receptor-binding domain (RBD) was performed
543 as previously described²⁶. The data were analyzed using a custom workflow implemented in
544 Snakemake⁶⁵. Briefly, paired-end reads were interleaved and merged into synthetic reads spanning
545 the entire RBD PCR amplicon using bbmerge.sh (v38.93) from the bbtools package (sourceforge.
546 net/projects/bbmap/) with default parameters. The merged reads were mapped to the SARS-CoV-2
547 reference sequence (Genbank MN908947.3) using minimap2 (v2.24) with the ``-ax sr`` preset for short
548 reads. The resulting mapping file was sorted with samtools (v1.14). Reads that fully contain the desired
549 amplicon sequence were extracted with the bedtools (v2.30.0) intersect tool. These reads were then
550 downsampled to a target depth of 1000 reads using reformat.sh (v38.93) from the bbtools package.
551 These downsampled reads were remapped to the MN908947.3 reference with minimap2. Residual
552 PCR primer sequences were then trimmed with samtools ampliconclip using the ``--hard-clip --both-
553 ends`` parameters. Next, a consensus sequence was generated by first generating a pileup with the
554 samtools mpileup tool using default settings and then generating a consensus with ivar (v1.3.1) using
555 the parameters ``-q 20 -t 0 -m 20``.

556 At the same time a consensus sequence was generated, the primer-trimmed reads were deduplicated
557 to determine how many of the reads were identical, essentially defining pseudo-haplotypes. Vsearch
558 (v2.21.1) fastx_uniques tool was used for deduplicating and enumerating the number of identical reads
559 in each sample.

560 Lineage-defining mutations in the RBD were used to differentiate Delta from Omicron consensus
561 sequences. Only one sample had evidence of mixed Delta and Omicron sequences.

562

563 *City of Milwaukee Health Department*

564 *Nucleic Acid Extraction and RT-qPCR*

565 Upon receipt of Thermo Scientific AerosolSense 2900 air sampler cartridge at City of Milwaukee
566 Health Department Laboratory (MHDL), collection substrates were aseptically removed and transferred
567 to a 5 mL sterile screw-cap tube filled with 1 mL of Remel viral transport medium. Samples were
568 kept frozen at -70°C until total nucleic extraction was performed using 200 µL elute and Applied
569 Biosystems™ MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit using ThermoFisher Scientific

570 KingFisher Flex instrument. Real-time RT-PCR setup was performed using 10 μ L of extract and
571 approved Applied Biosystems TaqPath™ COVID-19 Combo Kit containing three primer/probe sets
572 specific to different SARS-CoV-2 genomic regions (open reading frame 1ab (ORF1ab), spike (S)
573 protein, and nucleocapsid (N) protein-encoding genes) and primers/probes for bacteriophage MS2
574 which served as internal process control for nucleic acid extraction. RT-PCR assay was performed
575 on Applied Biosystems 7500 Fast Dx Real-Time PCR System according to the TaqPath™ COVID-19
576 Combo Kit protocol. PCR results were interpreted using the Applied Biosystems COVID-19 Interpretive
577 Software, and results were reported as positive, negative, or inconclusive for the detection of SARS-
578 CoV-2 RNA following successful quality control checks.

579

580 *SARS-CoV-2 Sequencing*

581 Samples with Ct values below 30 were sequenced using the ARTIC protocol and the Illumina DNA
582 Prep library kit on a MiSeq instrument ([https://www.protocols.io/view/sars-cov-2-sequencing-on-](https://www.protocols.io/view/sars-cov-2-sequencing-on-illumina-miseq-using-arti-bssjneqn)
583 [illumina-miseq-using-arti-bssjneqn](https://www.protocols.io/view/sars-cov-2-sequencing-on-illumina-miseq-using-arti-bssjneqn)). Data generated using the Integrated DNA Technologies ARCTIC
584 V4 primer panel were analyzed using the Illumina® DRAGEN COVID Lineage App, which uses a
585 customized version of the DRAGEN DNA pipeline to perform Kmer-based detection of SARS-CoV-2.
586 The app aligns reads to a reference genome, calls variants, and generates a consensus genome
587 sequence. Lineage/clade assignments were also confirmed using NextClade ([https://clades.nextstrain.](https://clades.nextstrain.org/)
588 [org/](https://clades.nextstrain.org/), version 1.14.0) and Pangolin COVID-19 Lineage Assigner (<https://pangolin.cog-uk.io/>, version
589 3.1.20) by uploading obtained FASTA files^{66,67}. Consensus sequences generated and related metadata
590 for environmental samples were shared publicly on Global Initiative on Sharing All Influenza Data
591 (GISAID) (www.gisaid.org), the principal repository for SARS-CoV-2 genetic information.

592

593 *Mayo Clinic*

594 *Nucleic Acid Extraction and RT-qPCR*

595 Upon removal of the screw cap from the air sample cartridge in the biosafety level 2 cabinet, the
596 air cartridge substrate was removed with a pair of disposable sterile forceps and transferred into a
597 sterile tube containing 1 mL of phosphate-buffered saline (PBS). The tube was vortexed for 10 sec,
598 and 200 μ L of the sample was used for nucleic acid extraction and purification on the KingFisher
599 Flex magnetic particle processor (Thermo Fisher Scientific, Inc.) using the MagMAX Viral/Pathogen II

600 Nucleic Acid Isolation Kit (Thermo Fisher Scientific, Inc.) and MVP_Flex_200ul software program, with
601 an elution volume of 50 uL. SARS-CoV-2 sequence targets (ORF1ab, N, and S gene sequences) were
602 amplified and detected with the FDA-authorized TaqPath COVID-19 Combo Kit (Life Technologies
603 Corp., Pleasanton, CA) on the Applied Biosystems 7500 Fast Dx Real-Time PCR System (Life
604 Technologies Corp.) per assay manufacturer's instructions for use.

605

606 **Ethics statement**

607 The Institutional Review Board of the University of Wisconsin-Madison Health Sciences waived ethical approval of this work.

608 References

- 609 1. Lee, J. S. et al. Analysis of the initial lot of the CDC 2019–Novel Coronavirus (2019-nCoV) real-time
610 RT-PCR diagnostic panel. *PLoS One* **16**, e0260487 (2021).
- 611 2. Wu, K. J. ‘It’s Like Groundhog Day’: Coronavirus Testing Labs Again Lack Key Supplies. *The New*
612 *York Times* (2020).
- 613 3. Supply shortages impacting COVID-19 and non-COVID testing. *ASM.org* [https://asm.org/](https://asm.org/Articles/2020/September/Clinical-Microbiology-Supply-Shortage-Collecti-1)
614 *Articles/2020/September/Clinical-Microbiology-Supply-Shortage-Collecti-1* (2021).
- 615 4. Vogels, C. B. F. et al. SalivaDirect: A simplified and flexible platform to enhance SARS-CoV-2
616 testing capacity. *Med (N Y)* **2**, 263–280.e6 (2021).
- 617 5. CDC. COVID data tracker. Centers for Disease Control and Prevention [https://](https://covid.cdc.gov/covid-data-tracker/?ACSTrackingID=FCP_1_USCDC_2145-DM75072&ACSTrackingLabel=%5BProof%201%5D%2002.04.2022%20-%20COVID-19%20Data%20Tracker%20Weekly%20Review&deliveryName=FCP_1_USCDC_2145-DM75072)
618 [covid.cdc.gov/covid-data-tracker/?ACSTrackingID=FCP_1_USCDC_2145-](https://covid.cdc.gov/covid-data-tracker/?ACSTrackingID=FCP_1_USCDC_2145-DM75072&ACSTrackingLabel=%5BProof%201%5D%2002.04.2022%20-%20COVID-19%20Data%20Tracker%20Weekly%20Review&deliveryName=FCP_1_USCDC_2145-DM75072)
619 [DM75072&ACSTrackingLabel=%5BProof%201%5D%2002.04.2022%20-%20COVID-19%20](https://covid.cdc.gov/covid-data-tracker/?ACSTrackingID=FCP_1_USCDC_2145-DM75072&ACSTrackingLabel=%5BProof%201%5D%2002.04.2022%20-%20COVID-19%20Data%20Tracker%20Weekly%20Review&deliveryName=FCP_1_USCDC_2145-DM75072)
620 [Data%20Tracker%20Weekly%20Review&deliveryName=FCP_1_USCDC_2145-DM75072](https://covid.cdc.gov/covid-data-tracker/?ACSTrackingID=FCP_1_USCDC_2145-DM75072&ACSTrackingLabel=%5BProof%201%5D%2002.04.2022%20-%20COVID-19%20Data%20Tracker%20Weekly%20Review&deliveryName=FCP_1_USCDC_2145-DM75072) (2020).
- 621 6. Da Silva, C., Murphy, J. & Lovelace, B., Jr. U.S. reports record 1.3M Covid cases in a day as
622 hospitalizations soar. *NBC News* [https://www.nbcnews.com/news/us-news/us-reports-record-](https://www.nbcnews.com/news/us-news/us-reports-record-13m-covid-cases-day-hospitalizations-soar-rcna11736)
623 [13m-covid-cases-day-hospitalizations-soar-rcna11736](https://www.nbcnews.com/news/us-news/us-reports-record-13m-covid-cases-day-hospitalizations-soar-rcna11736) (2022).
- 624 7. Kalish, H. et al. Undiagnosed SARS-CoV-2 seropositivity during the first 6 months of the
625 COVID-19 pandemic in the United States. *Sci. Transl. Med.* **13**, (2021).
- 626 8. Fact sheet: The Biden Administration to begin distributing at-home, rapid COVID-19 tests to
627 Americans for Free. *The White House* [https://www.whitehouse.gov/briefing-room/statements-](https://www.whitehouse.gov/briefing-room/statements-releases/2022/01/14/fact-sheet-the-biden-administration-to-begin-distributing-at-home-rapid-covid-19-tests-to-americans-for-free/)
628 [releases/2022/01/14/fact-sheet-the-biden-administration-to-begin-distributing-at-home-rapid-](https://www.whitehouse.gov/briefing-room/statements-releases/2022/01/14/fact-sheet-the-biden-administration-to-begin-distributing-at-home-rapid-covid-19-tests-to-americans-for-free/)
629 [covid-19-tests-to-americans-for-free/](https://www.whitehouse.gov/briefing-room/statements-releases/2022/01/14/fact-sheet-the-biden-administration-to-begin-distributing-at-home-rapid-covid-19-tests-to-americans-for-free/) (2022).
- 630 9. Rubin, R. First It Was Masks; Now Some Refuse Testing for SARS-CoV-2. *JAMA* **324**, 2015–2016
631 (2020).
- 632 10. News Division. Biden administration to invest more than \$12 billion to expand COVID-19 testing.
633 *HHS.gov* [https://www.hhs.gov/about/news/2021/03/17/biden-administration-invest-more-than-](https://www.hhs.gov/about/news/2021/03/17/biden-administration-invest-more-than-12-billion-expand-covid-19-testing.html)
634 [12-billion-expand-covid-19-testing.html](https://www.hhs.gov/about/news/2021/03/17/biden-administration-invest-more-than-12-billion-expand-covid-19-testing.html) (2021).

- 635 11. Larremore, D. B. et al. Test sensitivity is secondary to frequency and turnaround time for
636 COVID-19 screening. *Sci Adv* **7**, (2021).
- 637 12. Paltiel, A. D., Zheng, A. & Walensky, R. P. Assessment of SARS-CoV-2 Screening Strategies
638 to Permit the Safe Reopening of College Campuses in the United States. *JAMA Netw Open* **3**,
639 e2016818 (2020).
- 640 13. Chin, E. T. et al. Frequency of Routine Testing for Coronavirus Disease 2019 (COVID-19) in High-
641 risk Healthcare Environments to Reduce Outbreaks. *Clin. Infect. Dis.* **73**, e3127–e3129 (2021).
- 642 14. Mina, M. J., Parker, R. & Larremore, D. B. Rethinking covid-19 test sensitivity - A strategy for
643 containment. *N. Engl. J. Med.* **383**, e120 (2020).
- 644 15. Newman, C. M. et al. Initial Evaluation of a Mobile SARS-CoV-2 RT-LAMP Testing Strategy. *J.*
645 *Biomol. Tech.* **32**, 137–147 (2021).
- 646 16. Denny, T. N. et al. Implementation of a Pooled Surveillance Testing Program for Asymptomatic
647 SARS-CoV-2 Infections on a College Campus - Duke University, Durham, North Carolina, August
648 2-October 11, 2020. *MMWR Morb. Mortal. Wkly. Rep.* **69**, 1743–1747 (2020).
- 649 17. Mukherjee, U. K. et al. Author Correction: Evaluation of reopening strategies for educational
650 institutions during COVID-19 through agent based simulation. *Sci. Rep.* **11**, 20667 (2021).
- 651 18. Holden, T. M. et al. Geographic and demographic heterogeneity of SARS-CoV-2 diagnostic testing
652 in Illinois, USA, March to December 2020. *BMC Public Health* **21**, 1105 (2021).
- 653 19. Grigsby-Toussaint, D. S., Shin, J. C. & Jones, A. Disparities in the distribution of COVID-19 testing
654 sites in black and Latino areas in new York City. *Prev. Med.* **147**, 106463 (2021).
- 655 20. Cordes, J. & Castro, M. C. Spatial analysis of COVID-19 clusters and contextual factors in New
656 York City. *Spat. Spatiotemporal Epidemiol.* **34**, 100355 (2020).
- 657 21. Tao, R., Downs, J., Beckie, T. M., Chen, Y. & McNelley, W. Examining spatial accessibility to
658 COVID-19 testing sites in Florida. *Ann. GIS* **26**, 319–327 (2020).
- 659 22. Credit, K. Neighbourhood inequity: Exploring the factors underlying racial and ethnic disparities
660 in COVID-19 testing and infection rates using ZIP code data in Chicago and New York. *Reg. sci.*
661 *policy pract.* **12**, 1249–1271 (2020).

- 662 23. Rader, B. et al. Geographic access to United States SARS-CoV-2 testing sites highlights
663 healthcare disparities and may bias transmission estimates. *J. Travel Med.* **27**, (2020).
- 664 24. CDC. National Wastewater Surveillance System (NWSS). Centers for Disease Control and
665 Prevention [https://www.cdc.gov/healthywater/surveillance/wastewater-surveillance/wastewater-](https://www.cdc.gov/healthywater/surveillance/wastewater-surveillance/wastewater-surveillance.html)
666 [surveillance.html](https://www.cdc.gov/healthywater/surveillance/wastewater-surveillance/wastewater-surveillance.html) (2022).
- 667 25. Harris-Lovett, S. et al. Wastewater Surveillance for SARS-CoV-2 on College Campuses: Initial
668 Efforts, Lessons Learned, and Research Needs. *Int. J. Environ. Res. Public Health* **18**, (2021).
- 669 26. Gregory, D. A., Wieberg, C. G., Wenzel, J., Lin, C.-H. & Johnson, M. C. Monitoring SARS-CoV-2
670 Populations in Wastewater by Amplicon Sequencing and Using the Novel Program SAM Refiner.
671 *Viruses* **13**, (2021).
- 672 27. Peccia, J. et al. Measurement of SARS-CoV-2 RNA in wastewater tracks community infection
673 dynamics. *Nat. Biotechnol.* **38**, 1164–1167 (2020).
- 674 28. Ahmed, W. et al. First confirmed detection of SARS-CoV-2 in untreated wastewater in Australia:
675 A proof of concept for the wastewater surveillance of COVID-19 in the community. *Sci. Total*
676 *Environ.* **728**, 138764 (2020).
- 677 29. San Diego wastewater surveillance. SEARCH [https://searchcovid.info/dashboards/wastewater-](https://searchcovid.info/dashboards/wastewater-surveillance/)
678 [surveillance/](https://searchcovid.info/dashboards/wastewater-surveillance/) (2021).
- 679 30. Xu, Y. et al. Characteristics of pediatric SARS-CoV-2 infection and potential evidence for
680 persistent fecal viral shedding. *Nat. Med.* **26**, 502–505 (2020).
- 681 31. Chen, Y. et al. The presence of SARS-CoV-2 RNA in the feces of COVID-19 patients. *J. Med. Virol.*
682 **92**, 833–840 (2020).
- 683 32. Zheng, S. et al. Viral load dynamics and disease severity in patients infected with SARS-CoV-2
684 in Zhejiang province, China, January–March 2020: retrospective cohort study. *BMJ* **369**, m1443
685 (2020).
- 686 33. Yang, W., Elankumaran, S. & Marr, L. C. Concentrations and size distributions of airborne influenza
687 A viruses measured indoors at a health centre, a day-care centre and on aeroplanes. *J. R. Soc.*
688 *Interface* **8**, 1176–1184 (2011).

- 689 34. Prussin, A. J., 2nd & Marr, L. C. Sources of airborne microorganisms in the built environment.
690 Microbiome **3**, 78 (2015).
- 691 35. Wang, C. C. et al. Airborne transmission of respiratory viruses. Science **373**, (2021).
- 692 36. Myatt, T. A. et al. Detection of airborne rhinovirus and its relation to outdoor air supply in office
693 environments. Am. J. Respir. Crit. Care Med. **169**, 1187–1190 (2004).
- 694 37. Institute of Medicine (US), National Research Council (US) Committee on Effectiveness of National
695 Biosurveillance Systems: Biowatch & the Public Health System. The BioWatch System. (National
696 Academies Press (US), 2011).
- 697 38. Horve, P. F., Dietz, L., Northcutt, D., Stenson, J. & Van Den Wymelenberg, K. Evaluation of a
698 bioaerosol sampler for indoor environmental surveillance of Severe Acute Respiratory Syndrome
699 Coronavirus 2. PLoS One **16**, e0257689 (2021).
- 700 39. Dietz, L. et al. Exploring Integrated Environmental Viral Surveillance of Indoor Environments: A
701 comparison of surface and bioaerosol environmental sampling in hospital rooms with COVID-19
702 patients. bioRxiv (2021) doi:10.1101/2021.03.26.21254416.
- 703 40. Parhizkar, H. et al. Quantifying human and environmental viral load relationships amidst mitigation
704 strategies in a controlled chamber with participants having COVID-19. Research Square (2021)
705 doi:10.21203/rs.3.rs-940891/v1.
- 706 41. Horve, P. et al. Longitudinal analysis of built environment and aerosol contamination associated
707 with isolated COVID-19 positive individuals. Research Square (2021) doi:10.21203/rs.3.rs-861942/
708 v1.
- 709 42. Lednicky, J. A. et al. Collection of SARS-CoV-2 Virus from the Air of a Clinic Within a University
710 Student Health Care Center and Analyses of the Viral Genomic Sequence. Aerosol Air Qual. Res.
711 **20**, 1167–1171 (2020).
- 712 43. Santarpia, J. L. et al. Aerosol and surface contamination of SARS-CoV-2 observed in quarantine
713 and isolation care. Sci. Rep. **10**, 12732 (2020).
- 714 44. TrueMark respiratory panel 2.0 - US.

- 715 45. Kahn, J. S. et al. Seroepidemiology of human bocavirus defined using recombinant virus-like
716 particles. *J. Infect. Dis.* **198**, 41–50 (2008).
- 717 46. Wisconsin DHS Health Alert #35: Increased Influenza Identified in Midwest Universities.
718 Wisconsin Department of Health Services [https://content.govdelivery.com/accounts/WIDHS/
719 bulletins/2fe8698](https://content.govdelivery.com/accounts/WIDHS/bulletins/2fe8698).
- 720 47. Smyth, D. S. et al. Tracking cryptic SARS-CoV-2 lineages detected in NYC wastewater. *Nat.*
721 *Commun.* **13**, 635 (2022).
- 722 48. COVID-19: Wisconsin Coronavirus wastewater monitoring network. Wisconsin Department of
723 Health Services <https://www.dhs.wisconsin.gov/covid-19/wastewater.htm> (2020).
- 724 49. Governor Evers, DHS announce more than \$175 million to support COVID-19 testing in schools.
725 Wisconsin Department of Health Services [https://www.dhs.wisconsin.gov/news/releases/042021.
726 htm](https://www.dhs.wisconsin.gov/news/releases/042021.htm) (2021).
- 727 50. Wisconsin Schools. <https://www.greatschools.org/wisconsin/>.
- 728 51. Wisconsin child care centers. <https://childcarecenter.us/state/wisconsin>.
- 729 52. CMS changes Medicare payment to support faster COVID-19 diagnostic testing. [https://www.
730 cms.gov/newsroom/press-releases/cms-changes-medicare-payment-support-faster-covid-19-
731 diagnostic-testing](https://www.cms.gov/newsroom/press-releases/cms-changes-medicare-payment-support-faster-covid-19-diagnostic-testing).
- 732 53. Classification of Omicron (B.1.1.529): SARS-CoV-2 variant of Concern. [https://www.who.int/news/
733 item/26-11-2021-classification-of-omicron-\(b.1.1.529\)-sars-cov-2-variant-of-concern](https://www.who.int/news/item/26-11-2021-classification-of-omicron-(b.1.1.529)-sars-cov-2-variant-of-concern).
- 734 54. Lednicky, J. A. et al. Viable SARS-CoV-2 in the air of a hospital room with COVID-19 patients. *Int.*
735 *J. Infect. Dis.* **100**, 476–482 (2020).
- 736 55. Fabian, P., McDevitt, J. J., Houseman, E. A. & Milton, D. K. Airborne influenza virus detection with
737 four aerosol samplers using molecular and infectivity assays: considerations for a new infectious
738 virus aerosol sampler. *Indoor Air* **19**, 433–441 (2009).
- 739 56. Pan, M., Lednicky, J. A. & Wu, C.-Y. Collection, particle sizing and detection of airborne viruses. *J.*
740 *Appl. Microbiol.* **127**, 1596–1611 (2019).


- 741 57. Haktanir, A., Can, N., Seki, T., Kurnaz, M. F. & Dilmaç, B. Do we experience pandemic fatigue?
742 current state, predictors, and prevention. *Curr. Psychol.* 1–12 (2021) doi:10.1007/s12144-021-
743 02397-w.
- 744 58. Riemersma, K. K. et al. Shedding of infectious SARS-CoV-2 despite vaccination. *bioRxiv* (2021)
745 doi:10.1101/2021.07.31.21261387.
- 746 59. Andrews, N. et al. Covid-19 Vaccine Effectiveness against the Omicron (B.1.1.529) Variant. *N.*
747 *Engl. J. Med.* (2022) doi:10.1056/NEJMoa2119451.
- 748 60. Galani, A. et al. SARS-CoV-2 wastewater surveillance data can predict hospitalizations and ICU
749 admissions. *Sci. Total Environ.* **804**, 150151 (2022).
- 750 61. AerosolSense™ Sampler. [https://www.thermofisher.com/order/catalog/product/
751 AEROSOLSENSE?SID=srch-srp-AEROSOLSENSE](https://www.thermofisher.com/order/catalog/product/AEROSOLSENSE?SID=srch-srp-AEROSOLSENSE).
- 752 62. Aptima SARS-CoV-2 assay - Letter of Authorization. [https://www.fda.gov/media/138097/
753 download](https://www.fda.gov/media/138097/download).
- 754 63. CDC 2019-Novel Coronavirus (2019-nCoV) Real-Time RT-PCR Diagnostic Panel. [https://www.fda.
755 gov/media/134922/download](https://www.fda.gov/media/134922/download).
- 756 64. 2019-Novel Coronavirus (2019-nCoV) Real-time rRT-PCR Panel Primers and Probes. [https://www.
757 cdc.gov/coronavirus/2019-ncov/downloads/rt-pcr-panel-primer-probes.pdf](https://www.cdc.gov/coronavirus/2019-ncov/downloads/rt-pcr-panel-primer-probes.pdf).
- 758 65. Mölder, F. et al. Sustainable data analysis with Snakemake. *F1000Res.* **10**, 33 (2021).
- 759 66. Aksamentov, I., Roemer, C., Hodcroft, E. & Neher, R. Nextclade: clade assignment, mutation
760 calling and quality control for viral genomes. *J. Open Source Softw.* **6**, 3773 (2021).
- 761 67. O’Toole, Á. et al. Assignment of epidemiological lineages in an emerging pandemic using the pan-
762 golin tool. *Virus Evol* **7**, veab064 (2021).
- 763
764
765
766
767

768 **Acknowledgements**


769 This work was made possible by financial support through Rockefeller Regional Accelerator for
770 Genomic Surveillance (133 AAJ4558), Wisconsin Department of Health Services Epidemiology and
771 Laboratory Capacity funds (144 AAJ8216), and National Institutes of Health grant (U01DA053893). We
772 would like to thank all of the participating congregate settings for their partnership during this study.
773 We would like to thank Thermo Fisher Scientific, Inc. for donating 4 of the AerosolSense samplers
774 used in the data described in this work. We would like to acknowledge Elizabeth Durkes and Hannah
775 Kraussel, Public Health Emergency Response Planning Coordinators, Lindor Schmidt, Environmental &
776 Disease Control Specialist, and Dr. Nicholas Tomaro, Emergency Preparedness Environmental Health
777 Director at the City of Milwaukee Health Department for their assistance in coordinating site selection,
778 sampling, cartridge replacement, and unit decontamination processes.

778 **Contributions**

780 M.D.R contributed to the conceptualization, data curation, formal analysis, investigation, methodology,
781 project administration, visualization, writing—original draft preparation, writing—review and editing.
782 C.M.N. contributed to data curation, investigation, methodology, supervision, writing—review
783 and editing. S.F.B. contributed to project administration, data curation, resources, visualization,
784 writing—review and editing. M.R.S. contributed to the investigation, data curation, methodology,
785 writing—review and editing. R.W.W. contributed to data curation, methodology, writing—review
786 and editing. A.K.Y. contributed to resources, writing—review and editing. E.J.O. contributed to
787 software, data curation, writing—review and editing. N.D. contributed to the investigation, data
788 curation, methodology, writing—review and editing. A.L. contributed to the investigation, formal
789 analysis, data curation, methodology, supervision, writing—review and editing. K.P.P. contributed to
790 resources, supervision, writing—review and editing. N.S. contributed to resources, supervision, data
791 curation, writing—review and editing. J.A.M. contributed to resources, supervision, data curation,
792 writing—review and editing. M.A.A. contributed to resources, supervision, data curation, writing—
793 review and editing. W.M.R contributed to resources, supervision, data curation, writing—review and
794 editing. J.A.Z. Contributed to the investigation, developed methods, sample processing, performed
795 analysis, data curation, writing—review and editing. M.K. Contributed to the method development,
796 analysis, data curation, writing—review and editing. L.J.B. contributed to conceptualization, data
797 curation, formal analysis, investigation, methodology, writing—review and editing. E.C.B. contributed
798 to the investigation, formal analysis, data curation, methodology, writing—review and editing. D.C.
799 contributed to the investigation, formal analysis, data curation, methodology, writing—review and
800 editing. C.R. contributed to the investigation, data curation, methodology, writing—review and editing.
801 D.A.G contributed to formal analysis, data curation, writing—review and editing. J.D.Y. contributed
802 to conceptualization, supervision, methodology, writing—review and editing. S.B. contributed
803 to resources, conceptualization, supervision, writing—review and editing. M.C.J. contributed to
804 the investigation, data curation, methodology, writing—review and editing. M.T.A. contributed to
805 conceptualization, supervision, methodology, writing—review and editing. T.C.F. contributed to
806 the conceptualization, funding acquisition, writing—review and editing. D.H.O. contributed to the
807 conceptualization, formal analysis, software, funding acquisition, methodology, supervision, project
808 administration, writing—original draft preparation, writing—review and editing. S.L.O. contributed to
809 the conceptualization, funding acquisition, methodology, supervision, project administration, writing—
810 review and editing.




811 **Figure 01. Air sample testing workflow.** (A) Overview of air sample collection, processing, and testing. (B) Air sample data
 812 collection and management. Individuals in charge of changing air cartridges at surveillance sites use the iOS and Android
 813 Askidd mobile app to collect metadata on air samples when cartridges are inserted and removed. Data are compiled in
 814 Labkey database and displayed to surveillance sites in the Askidd mobile app. Created with BioRender.com.




815 **Figure 02. COVID-19 outbreak timeline.** Confirmed COVID-19 cases and air sample SARS-CoV-2 RT-qPCR results in the
816 congregate setting are plotted over time. Orange dots represent confirmed COVID-19 cases from individuals present in the
817 building. Blue boxes show the number of COVID-19 cases that occurred while close contacts were in quarantine. Air sample
818 SARS-CoV-2 RT-qPCR results are represented by boxes as positive (green), negative (red), or inconclusive (orange). The gray
819 box indicates that no sample was collected during that time period.

Sampler/Day	Day 1	Day 2	Day 3	Day 4	Day 5	Day 6	Day 7
Continuous Wk 1	35.4 36.5						
Daily Wk 1	Neg	37.4 37.1	Neg	Neg	Neg	Neg	Neg
Continuous Wk 2	36.4 Undet.						
Daily Wk 2	36.6 Undet.		Neg	Neg	Neg	Neg	Neg
Continuous Wk 3	Neg						
Daily Wk 3	Neg	Neg			Neg	Neg	Neg
Continuous Wk 4	Neg						
Daily Wk 4	Undet. 39.1	Neg	Neg	Neg	Neg	Neg	Neg
Continuous Wk 5	Neg						
Daily Wk 5	Neg	Neg	Neg	Neg	Neg	Neg	Undet. 39.4

820 **Figure 03. Comparison of SARS-CoV-2 RT-qPCR results from continuous and daily air sampling intervals.** Two adjacent
821 Thermo Scientific AerosolSense instruments were run continuously or daily over several days. SARS-CoV-2 genomic material
822 was detected by two separate RT-qPCR CDC assays. If the cycle threshold values of one or both N1 and N2 assays were
823 less than 40, the Ct values are shown in the box separated by a '|'. 'Undet.' was used for assays that had Ct values greater
824 than 40. If the Ct values of both RT-qPCR assays were greater than 40, the boxes are labeled as 'Neg'. Samples considered
825 to be positive are shaded red. Boxes shaded in grey are either inconclusive or negative.



826 **Figure 04. In-air respiratory pathogen detection in congregate settings.** (A) Respiratory pathogen detection in air samples
 827 collected from a K-12 school, (B) preschool, (C) campus coffee shop, and (D) campus athletic facility. Genomic material from
 828 40 respiratory pathogens was detected by semi-quantitative RT-PCR using the TrueMark Respiratory 2.0 TaqMan Array Card.
 829 SARS-CoV-2 genomic material was detected by two separate RT-qPCR CDC assays. Boxes shaded in red, pink, and gray
 830 represent positive, inconclusive, and negative air samples collected during the sampling interval on the x-axis. No sample
 831 was tested for boxes shaded in black.



832 **Figure 05. Detection of SARS-CoV-2 and influenza A virus in Dane County, WI.** (A) Influenza A virus (IAV) detection in
 833 air samples collected from congregate settings. IAV genomic material was detected by semi-quantitative RT-PCR using the
 834 TrueMark Respiratory 2.0 TaqMan Array Card. (B) SARS-CoV-2 detection in air samples collected from congregate settings.
 835 SARS-CoV-2 genomic material was detected by two separate RT-qPCR N1 and N2 CDC assays. Boxes shaded in red, pink,
 836 and gray represent positive, inconclusive, and negative air samples collected during the sampling interval in the x-axis.
 837 No sample was tested for boxes shaded in black. Campus sites were located on the college campus of the University of
 838 Wisconsin-Madison.

Table 1. SARS-CoV-2 air sample results.

Location	Site Name	Start Date	End Date	Number of Samples	Positive	Negative	Inconclusive
Dane County, WI	Preschool #1	8/18/21	2/8/22	49	3	43	3
	Preschool #2	8/11/21	10/14/21	22	2	18	2
	K-12 School #1	7/26/21	2/8/22	73	4	62	7
	K-12 School #2	10/14/21	2/9/22	15	8	5	2
	K-12 School #3	12/14/21	2/8/22	7	7	0	0
	K-12 School #4	12/15/21	2/8/22	8	6	1	1
	Hospital	8/20/21	10/25/21	51	18	33	0
	Campus Athletic Facility	7/19/21	2/9/22	179	20	141	18
	Campus Coffee Shop	8/17/21	2/3/22	54	5	44	5
	Office	9/30/21	12/10/21	8	0	8	0
Minneapolis, MN	Brewery taproom	10/18/21	2/7/22	26	11	2	13
Rochester, MN	Bar	9/27/21	11/24/21	9	5	4	0
	Hospital Cafeteria	9/20/21	11/24/21	10	6	4	0
Milwaukee, WI	Emergency Housing Facility #1	12/17/21	2/8/22	9	5	3	1
	Emergency Housing Facility #2	12/17/21	2/8/22	7	6	1	0
Total				527	106	369	52

Table 2. Air sample SARS-CoV-2 sequencing results.

Location	Air sample barcode	Start	Finish	SARS-CoV-2 N gene PCR Ct	Lineage(s)	Spike RBD amino acid differences vs. SARS-CoV-2 reference	Accession Number
Emergency Housing Facility #1	AE000010795F4A	12/21/21	1/7/22	31.77	BA.1	S477N; T478K; E484A; Q493R; G496S; Q498R; N501Y; Y505H; T547K	EPI_ISL_8879389
Emergency Housing Facility #2	AE000010795B42	12/21/21	1/7/22	25.94	BA.1	S477N; T478K; E484A; Q493R; G496S; Q498R; N501Y; Y505H; T547K	EPI_ISL_8879388
Brewery taproom	AE000010466C36	11/22/21	11/29/21	30.5 34.28	Delta	L452R; T478K	SRX14331279
Brewery taproom	AE000010464938	11/22/21	11/29/21	31.63 35.15	Delta	L452R; T478K	SRX14331280
Brewery taproom	AE000010467837	12/6/21	12/13/21	38 69.63	Delta	L452R; T478K	SRX14331281
Brewery taproom	AE000010467A3C	12/6/21	12/13/21	37.41 42.9	Delta	L452R; F456L; T478K; F562F	SRX14331282
Brewery taproom	AE000010463B32	12/30/21	1/3/22	35.28 38.77	BA.1	K417N; N440K; G446S; S477N; T478K; E484A; Q493R; G496S; Q498R; N501Y; Y505H; T547K	SRX14331283
Brewery taproom	AE00001046442E	12/30/21	1/3/22	35.69 37.9	Delta and BA.1	K417N; N440K; G446S; L452R; S477N; T478K; E484A; Q493R; G496S; Q498R; N501Y; Y505H; T547K	SRX14331284
Brewery taproom	AE000010463F3A	1/3/22	1/10/22	33.85 37.47	BA.1	K417N; N440K; G446S; S477N; T478K; E484A; Q493R; G496S; Q498R; N501Y; Y505H; T547K	SRX14331285
Brewery taproom	AE000010465530	1/10/22	1/17/22	33.37 36.34	BA.1	K417N; N440K; G446S; S477N; T478K; E484A; Q493R; G496S; Q498R; N501Y; Y505H; T547K	SRX14331286
Brewery taproom	AE00001053FA3D	1/17/22	1/25/22	34.04 38.94	BA.1	K417N; N440K; G446S; S477N; T478K; E484A; Q493R; G496S; Q498R; N501Y; Y505H; T547K	SRX14331287

Samples with two cycle threshold (Ct) values listed in the table were tested with two SARS-CoV-2 N1 and N2 assays. N1 and N2 Ct values are separated by '|'. Samples with one Ct value listed in the table were tested with the Applied Biosystems TaqPath™ COVID-19 Combo Kit. Abbreviations: RBD, receptor binding domain; Ct, cycle threshold.