

SUPPLEMENTAL MATERIAL

This appendix has been provided by the authors to give readers additional information about their work.

Title: Cardiovascular Signatures of COVID-19 Predict Mortality and Identify Barrier Stabilizing Therapies

AUTHOR NAMES: Dakota Gustafson, B.Sc.H.^{a,b}, Michelle Ngai, Ph.D.^a, Ruilin Wu, M.Sc.^{a,b}, Huayun Hou, Ph.D.^c, Alice Schoffel, M.P.H.^d, Clara Erice, Ph.D.^e, Serena Mandla, M.Sc.^f, Filio Billia, M.D., Ph.D.^{a,d}, Michael D. Wilson, Ph.D.^{c,g}, Milica Radisic, Ph.D.^f, Eddy Fan, M.D., Ph.D.^{a,h,i}, Uriel Trahtemberg, M.D., Ph.D.^{j,k}, Andrew Baker, M.D.^{h,i,k}, Chris McIntosh, Ph.D.^{d,l,m,n,o}, Chun-Po S. Fan, Ph.D.^d, Claudia C. dos Santos, M.D.^{h,j}, Kevin C. Kain, M.D.^{a,b}, Kate Hanneman, M.D., MPH^{a,d,l}, Paaladinesh Thavendiranathan, M.D., S.M.^{d,i,l,p}, Jason E. Fish, Ph.D.^{a,b,d,i,*}, Kathryn L. Howe, M.D., Ph.D.^{a,d,i,q*}

AFFILIATIONS: ^aToronto General Hospital Research Institute, University Health Network, Toronto, Canada; ^bDepartment of Laboratory Medicine and Pathobiology, University of Toronto, Toronto, Canada; ^cProgram in Genetics and Genome Biology, The Hospital for Sick Children, Toronto, Canada; ^dPeter Munk Cardiac Centre, Toronto General Hospital, University Health Network, Toronto, Canada; ^eJohns Hopkins School of Medicine, Baltimore, USA; ^fInstitute of Biomedical Engineering, University of Toronto, Toronto, Canada; ^gDepartment of Molecular Genetics, University of Toronto, Toronto, Canada; ^hInterdepartmental Division of Critical Care and Institute of Medical Sciences, University of Toronto, Toronto, Canada; ⁱInstitute of Medical Science, University of Toronto, Toronto, Canada; ^jKeenan Research Center for Biomedical Research, Unity Health Toronto, Toronto, Canada; ^kCritical Care Department, Galilee Medical Center, Nahariya, Israel; ^lJoint Department of Medical Imaging, University Health Network, University of Toronto, Toronto, Canada; ^mTechna Institute, University Health Network, Toronto, Canada; ⁿDepartment of Medical Biophysics, University of Toronto, Toronto, Canada; ^oVector Institute, University of Toronto, Toronto, Canada; ^pTed Rogers Program in Cardiotoxicity Prevention, Toronto General Hospital, Toronto, Canada; ^qDivision of Vascular Surgery, Department of Surgery, University of Toronto, Toronto, Canada

TABLE OF CONTENTS

SUPPLEMENTAL ABBREVIATIONS AND ACRONYMS.....	4
SUPPLEMENTAL METHODS.....	6
Rationale and Design.....	6
Patient Categorization.....	6
Study Population.....	6
Clinical Data Collection.....	6
Processing of Bloodwork.....	7
Protein Biomarker Analysis.....	7
Cell Culture.....	7
xCelligence Real-Time Cell Analysis (Transendothelial Electrical Resistance [TEER]).....	7
Transwell Leak Assay.....	8
Bulk RNA-Seq.....	8
HTG EdgeSeq MicroRNA (miRNA) Whole Transcriptome Assay (WTA) from Plasma.....	9
HTG EdgeSeq MiRNA Analysis.....	9
Simultaneous Multiplexed <i>In Vitro</i> Immunofluorescence.....	9
SUPPLEMENTAL TABLES.....	10
Online Table I. Key Reagents and Resources.....	10
Online Table II. Antibodies.....	12
Online Table III. Admission Clinical Laboratory Findings.....	13
Online Table IV. Association of Baseline Clinical Characteristics to Mortality Amongst All Admitted COVID-19 Patients.....	14
Online Table V. Intensive Care Unit-Level Clinical Characteristics and Clinical Outcomes.....	16
Online Table VI. Association of Baseline Clinical Characteristics to Mortality Amongst All Severe COVID-19 Patients.....	17
SUPPLEMENTAL FIGURES.....	19
Online Figure I; Related to Methods and Table 1: Flow diagram of patients enrolled between the COLOBILI Study (St. Michael's Hospital) and the COVID Study (University Health Network).....	19
Online Figure II; Related to Methods and Table 1: Spike (trimer) antigen serology testing from patients having a negative SARS-CoV-2 polymerase chain reaction result.....	20
Online Figure III; Related to Figure 1: The association of coronary artery disease with mortality characterized in terms of proportion of deceased patients stratified by status.....	21
Online Figure IV; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers amongst the entire cohort (SARS-CoV-2 negative and positive populations).....	22
Online Figure V; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the mild COVID-19 subgroup.....	23
Online Figure VI; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the moderate COVID-19 subgroup.....	24
Online Figure VII; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the severe COVID-19 subgroup.....	25

Online Figure VIII; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the mild SARS-CoV-2 negative subgroup.....	26
Online Figure IX; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the severe SARS-CoV-2 negative subgroup.....	27
Online Figure X; Related to Figure 2. Plasma Concentration of Endothelial Dysfunction and Inflammatory Markers at t_{0-1}	28
Online Figure XI; Related to Figure 2: Plasma Concentration of Endothelial Dysfunction and Inflammatory Markers at t_{0-1} and ability to discriminate survival in ICU patients.	31
Online Figure XII; Related to Figure 2: Plasma Concentration of Endothelial Dysfunction and Immunological Markers at t_{0-1} in ICU patients.....	33
Online Figure XIII; Related to Figure 2: Plasma Concentration of (a) IL-6 and (b) MPO, longitudinally between severe COVID-19 patients and severe SARS-CoV-2 negative patients.....	34
Online Figure XIV; Related to Figure 3 and 4: Plasma MicroRNA Transcriptome Across the Disease Severity Subgroups.	36
Online Figure XV; Related to Figure 3 and 4: Plasma MicroRNA Transcriptome Across the Disease Severity Subgroups.	38
Figure XVI; Related to Figure 4: Feature importance of a machine learning model incorporating clinical data. ..	39
Online Figure XVII; Related to Figure 4: Feature importance of a machine learning model incorporating both clinical data and protein expression metrics.	40
Online Figure XVIII; Related to Figure 5: Association of Biomarkers with In-Hospital Mortality for Severe COVID-19 Patients.	41
Online Figure XIX; Related to Figure 5: T_{0-1} COVID-19 Patient Plasma Selectively Induces Acute Increases in Endothelial Permeability.....	42
Online Figure XX; Related to Figure 5: Correlation of t_{0-1} Plasma Cardiovascular Biomarkers in COVID-19 positive patients to Induction of Endothelial Permeability.....	43
Online Figure XXI; Related to Figure 6: Endogenous sSlit2 is upregulated in severe COVID-19 patient plasma.	44
SUPPLEMENTAL DATA FILE ANNOTATIONS.....	45
Supplementary Data File I. Quality control table for all RNA-sequencing experiments used in this study.	45
Supplementary Data File II. R documentation file for the analysis of the RNA-sequencing experiments.	45
Supplementary Data File III. R code for the analysis of the RNA-sequencing experiments.	45
Supplementary Data File IV. Full list of differentially expressed miRNA with pairwise comparisons between COVID-19 cohorts and the negative controls.....	45
Supplementary Data File V. Full list of pathway enrichments for miRNA-sequencing experiment between COVID-19 cohorts and the negative controls.....	45
Supplementary Data File VI. Full list of differentially expressed genes with pairwise comparisons between COVID-19 cohorts and the negative controls.....	45
Supplementary Data File VII. Full list of pathway enrichments for mRNA-sequencing experiment between COVID-19 cohorts and the negative controls.....	45
Supplementary Data File VIII. Gene set enrichment analysis for mRNA-sequencing experiment between COVID-19 cohorts and the negative controls.....	45
SUPPLEMENTAL REFERENCES	46

SUPPLEMENTAL ABBREVIATIONS AND ACRONYMS

ACE	Angiotensin-converting-enzyme inhibitors
Ang-2	Angiopoietin-2
APACHE	Acute physiologic assessment and chronic health evaluation
ARBs	Angiotensin II receptor blockers
ARDS	Acute respiratory distress syndrome
BMI	Body mass index
CAD	Coronary artery disease
Cat#	Catalogue number
CCB	Calcium channel blocker
CF	Cystic fibrosis
CKD	Chronic kidney disease
COPD	Chronic obstructive pulmonary disease
COVID-19	Coronavirus Disease 2019
CV	Cardiovascular
DAPI	4',6-diamidino-2-phenylindole
EBM	Endothelial basal media
ECMO	Extracorporeal membrane oxygenation
ED	Emergency department
EGM	Endothelial growth medium
ET-1	Endothelin-1
FDR	False discovery rate
FITC	Fluorescein isothiocyanate
FiO₂	Fraction of inspired oxygen
GERD	Gastroesophageal reflux disease
GSEA	Gene set enrichment analysis
GSE	Gene expression omnibus series accession number
Hs-cTnI	High-sensitivity cardiac troponin
Ig	Immunoglobulin
IL	Interleukin
IQR	Interquartile range
kDa	Kilodalton
KEGG	Kyoto encyclopedia of genes and genomes
miR/miRNA	MicroRNA
MIV	Mechanical ventilation
MPO	Myeloperoxidase
NIV	Non-invasive ventilation
NSTEMI	Non-ST-elevation myocardial infarction
OSA	Obstructive sleep apnea
PAD	Peripheral artery disease
PF	Ratio of arterial oxygen partial pressure to fractional inspired oxygen
PFO	Patent foramen ovale
pHUVEC	Pooled human umbilical vein endothelial cells
QC	Quality control
REB	Research ethics board
RR	Respiratory rate
RRID:AB	Research resource identifier (antibody)
RTCA	Real-time cell analyzer
S.D.	Standard deviation
SARS-COV-2	Severe acute respiratory syndrome coronavirus-2
sICAM-1	Soluble intercellular adhesion molecule-1
SOFA	Sequential organ failure assessment
sSLIT2	Soluble slit guidance ligand 2
sTREM-1	Soluble triggering receptor expressed on myeloid cells-1

sVCAM-1	Soluble vascular cell adhesion molecule-1
TIA	Transient ischemia attack
TC	Tissue culture
TEER	Transendothelial electrical resistance
TIA	Transient ischemic attack
TNFα	Tumor necrosis factor alpha
VTE	Venous thromboembolism
WBCs	White blood cells

SUPPLEMENTAL METHODS

Rationale and Design: Coronavirus Disease 2019 (COVID-19) is a leading infectious etiology currently present throughout the world, for which there are limited therapeutic interventions. Understanding the pathobiology of COVID-19 outcomes can enable personalized patient management protocols, improve survival, and aid in the preparation of diagnostic tools for future pandemics. Due to the dynamic nature of clinical care guidelines, clinical care did not follow a standardized protocol but rather was determined by individual providers and patient needs. Plasma of enrolled patients was taken at presentation and at set intervals (days 2-3 [t₂₋₃], 4-5 [t₄₋₅], and 6-7 [t₆₋₇]) following obtainment of consent. Comprehensive, integrated analysis of plasma transcriptome data was performed to prioritize COVID-19 outcome signals.

Patient Categorization: After retrospective clinical adjudications of medical records and triaging protocols, patients were assigned to pre-defined clinical groups, centered around the National Institutes of Health, 'Clinical Spectrum of SARS-CoV-2 Infection'(1), those being: (i) severe acute respiratory syndrome coronavirus (SARS-CoV-2) negative patients with mild acute respiratory disease, (ii) mild COVID-19, (iii) moderate COVID-19, (iv) severe COVID-19, and (v) SARS-CoV-2 negative patients requiring intensive care unit (ICU) level management for severe respiratory illness. These assignments were made solely on the basis of information available in the medical record and were blind to any novel biomarker data, which had not yet been generated. For biomarker studies, where appropriate, patients were matched for age, sex, body mass index (BMI), and co-existing conditions with non-infected controls as reference groups.

Study Population: This is a multicenter, secondary analysis of a prospectively recruited longitudinal cohort study enrolling consecutive patients with suspected SARS-CoV-2 infection who were referred to two Canadian quaternary care networks in Toronto, Canada from May 2020 to December 2020: University Health Network and St. Michael's Hospital. All participants who were 18 years of age or older, provided either direct written informed consent or were consented into the study by a lawfully entitled substitute decision-maker on behalf of a participant when lacking the capacity to make the decision. The study, and consenting, was conducted in accordance with protocols approved by the Research Ethics Board (REB) of the University Health Network (REB#: 20-5453.6; Cardiovascular Disease and Outcomes among Patients with SARS-CoV-2 Infection During Admission and Post-Discharge [The COVID study]) or St. Michael's Hospital (REB#: 20-078; COVID-19 Longitudinal Biomarkers in Lung Injury [COLOBILI]). SARS-CoV-2-negative patients with severe respiratory illness symptoms were enrolled within the COLOBILI study. Diagnoses of COVID-19 were confirmed through real-time reverse transcription-polymerase chain reaction assays of nasopharyngeal swabs according to the Public Health Ontario guidelines for SARS-CoV-2 testing.(2) The timeframe for recruitment largely excludes community level spread for the variants of concern (i.e., B.1.1.7 [Alpha], B.1.351 [Beta], B.1.617.2 [Delta], P.1 [Gamma], and B.1.1.529 [Omicron]), with the assumption that all patients harbored one of the predominant unmutated G strains (i.e., GR, GH, and GV).(3) Samples were collected prior to the initiation of public vaccination programs in Ontario, with all patients assumed to be unvaccinated. Admitted patients were followed up after their COVID-19 diagnosis, with all causes of in-hospital mortality, complications, and therapeutic regimen ascertained until discharge. Patients' data were extracted from the in-hospital electronic medical records, de-identified, and assigned random identification numbers which were used throughout the project. Information on sex, age, and pertinent clinical parameters for the cohorts are provided (Table 1). The laboratory parameters were collected as reported by the individual centers, with standard international reference ranges applied to decide the cutoff point for abnormal levels.

Clinical Data Collection: Clinical characteristics, medical history, therapeutics administered during admission, complications, and outcomes were obtained and extracted from electronic medical records by clinical coordinators. Participant disease severity was quantified according to the National Institutes of Health Clinical Spectrum of SARS-CoV-2 infection which is largely based on respiratory parameters. Complete follow-up was available only for those requiring admission to hospital. Obesity was defined as BMI ≥ 30 kg/m² in line with World Health Organization guidelines.(4) Preexisting cardiovascular disease (i.e., atrial fibrillation, arterial hypertension, coronary artery disease, dyslipidemia, diabetes mellitus, chronic obstructive pulmonary disease, or heart failure), myocardial infarction, and concurrent active malignancy were defined and recorded from the electronic medical records at the discretion of the clinical coordinators upon availability from previously recorded histories. Smoking was collected as a self-reported variable, with social smoking being defined as less than <4 times per week, and avid smoker being anything greater than that. Cardiovascular complications were defined as a surrogate of new-onset atrial or ventricular arrhythmia (atrial fibrillation, atrial flutter, non-sustained ventricular tachycardia, sustained ventricular tachycardia), acute

coronary syndrome, clinical heart failure, right ventricular failure, left ventricular failure, biventricular failure, moderate or greater pulmonary embolism, tamponade, stroke, acute non-coronary ischemia due to hypercoagulable state, or intracardiac thrombus. In patients admitted to the ICU the Acute Physiologic Assessment and Chronic Health Evaluation II and the Sequential Organ Failure Assessment scores were used.(5, 6) These are validated methods for grading the severity of illness in critically ill patients based upon point-scoring systems associated with the degree of dysfunction with specific physiologic variables.

Processing of Bloodwork: Peripheral blood samples were collected synchronously with standard of care bloodwork at the University Health Network or St. Michael's Hospital between May 2020 and December 2020. Peripheral blood samples (10 mL) were drawn from the cubital vein into BD Vacutainer® Blood Collection Tubes (BD Bioscience, Franklin Lakes, NJ) containing K₂EDTA and processed within three hours. Plasma was separated from whole blood through centrifugation (2,000×g, 24°C, 15 min) and stored at -80°C until downstream processing. At no time during the process was the plasma subjected to temperatures below 4°C or above 25°C. Samples were thawed on ice and the plasma was subjected to sequential centrifugation of (2,500×g, 4°C, 25 min) to reduce platelets and large particulate. Hemolysis was examined prior to downstream analysis by measuring the absorbance at 414 nm using a DS-11⁺ Spectrophotometer (DeNovix, Wilmington, Delaware, United States).

Protein Biomarker Analysis: Circulating levels of angiopoietin-2 (Ang-2; Lower limit of quantification [LLOQ] 9.91 pg/mL), soluble CD62 antigen-like family member E (sE-Selectin; LLOQ 4.22 pg/mL), soluble CD54/intercellular adhesion molecule-1 (sICAM-1; LLOQ 4.1 pg/mL), soluble CD106/vascular cell adhesion protein-1 (sVCAM-1; LLOQ 137 pg/mL), CD105/endoglin, endothelin-1 (ET-1; LLOQ 0.250 pg/mL), interleukin-6 ([IL]-6; LLOQ 0.41 pg/mL), IL-8 (LLOQ 0.19 pg/mL), and soluble triggering receptor expressed on myeloid cells-1 (sTREM-1; LLOQ 4.19 pg/mL) were quantified in platelet free plasma samples using the Simple Plex Ella (ProteinSimple, San Jose, CA, USA) multiplex platform according to the manufacturer's instructions; all Simple Plex values are reported as the average of triplicate readings. Soluble Slit homolog 2 protein (sSLIT2) was quantified using an enzyme-linked immunosorbent sandwich assay (ELISA; LLOQ 0.10 ng/mL, Elabscience, Wuhan, China). Myeloperoxidase was additionally quantified through ELISA (LLOQ; 0.062 ng/mL; R&D Systems Minneapolis, MN, USA). Circulating cardiac troponin I (cTnI) was quantified using a CLIA certified high sensitivity ELISA kit (LLOQ 0.92 pg/mL; Biomatik, Kitchener, ON, CAN). Cardiac injury was defined as plasma levels of high-sensitivity cTnI (hs-cTnI) greater than the 99th percentile of normal values, as per clinical guidelines. See Online Table I for reagents.

Cell Culture: Pooled human umbilical vein endothelial cells ([pHUECs], Lonza, Basel, Switzerland) from multiple donors were cultured in EC Growth Basal Medium-2 (Lonza, Basel, Switzerland) containing the complete EC Growth Medium BulletKit™ (Lonza, Basel, Switzerland) at 37°C in a 5% CO₂ humidified incubator. Cells were maintained on 10 cm tissue culture plates (Corning, Costar, New York, NY, USA) coated with 0.1% (v/v) gelatin attachment factor (ThermoFisher, Waltham, MA, USA) and passaged every 2-4 days at 70%-80% confluency. For experimentation, pHUECs at passages 3-7 were used. Cells tested negative for the presence of mycoplasma contamination (ThermoFisher, Waltham, MA, USA).

xCelligence Real-Time Cell Analysis (Transendothelial Electrical Resistance [TEER]): Baseline resistance measurements were acquired using 150 µL EC Growth Basal Medium-2 for 15 minutes prior to experiment initiation using xCelligence Real-Time Cell Analyzer (RTCA, ACEA Biosciences, San Diego, CA). Following this, pHUECs were seeded at densities of 4x10⁴/well into multiple 16-well E-plates (ACEA Biosciences, San Diego, CA) in volumes of 150 µL EC Growth Basal Medium-2 and monitored on the xCelligence RTCA. Once the cells reached confluency, 20% v/v plasma or 2 U/mL thrombin (Sigma-Aldrich, St Louis, MO) was added, and the cells incubated at 37°C in a 5% carbon dioxide atmosphere for six hours. To quantify the change in permeability, cell indexes (surrogate metric for change in resistance across the monolayer) were adjusted to the resistance of reference wells (i.e., untreated cells in growth media) and normalized to the time point immediately prior to the addition of the plasma. The net area between the curve of the normalized dataset was utilized to calculate the overall change in permeability, with an increased area representing increases in TEER and decreases representing declines in TEER (i.e., loss of endothelial barrier function). During therapeutic testing, cells were co-treated with 400 µM QHREDGS peptide (Q-Peptide)(7), 1280 ng/mL recombinant Slit2-N(8), 50 µg/mL nangibotide(9), 3 µM dexamethasone(10), or a matching DMSO control with Dulbecco's Phosphate Buffered Saline without magnesium and calcium ([PBS⁻], Gibco, Gaithersburg, MD, USA).

Transwell Leak Assay: Endothelial monolayer leak assays were performed as previously described.(11) Briefly, pHUVECs were seeded on 3 μ m pore transwell inserts (Corning Life Sciences, Corning, NY, USA) coated with 0.1% (v/v) gelatin attachment factor (ThermoFisher, Waltham, MA, USA). The cells were grown on the inserts for two days until reaching 90%-95% confluency. Cells were subsequently treated for either one hour (acutely) *or* six hours (sub-chronically) with either 20% (v/v) human plasma or control PBS^{-/-}. Prior to leak quantification, cell media was changed to Hanks Buffered Salt Solution (ThermoFisher, Waltham, MA, USA) and 1 mg/mL 40 kilodalton (kDa) fluorescein isothiocyanate–dextran (FITC-dextran, Sigma-Aldrich, St Louis, MO) was added to the top chamber with tracer flux subsequently being allowed to occur for one hour; 2 U/mL thrombin was added as a positive control to certain wells at the time of FITC addition. The experimenter was blinded to the grouping of each sample. The FITC accumulation in the bottom chamber was assessed in triplicate via excitation at 485 nm and emission at 535 nm using a Biotek Cytation 5 (Biotek, Winooski, VT, USA).

Bulk RNA-Seq: Total RNA was isolated from treated pHUVEC samples (n=5 control [PBS^{-/-}], n=5 mild SARS-CoV-2 negative, n=5 mild disease, n=5 moderate disease, and n=5 severe disease) at the six-hour timepoint using the RNeasy Plus Micro kit (Qiagen, Germantown, MD, USA), after washing twice with ice-cold PBS^{-/-}. RNA quantities and quality were assessed using the Agilent 2100 Bioanalyzer (Agilent Technologies, Mississauga, ON, Canada). All samples passed a quality control threshold (RNA integrity number ≥ 7.0) to proceed to library preparations and RNA-seq. A total amount of 20 ng RNA per sample was used as input material for the RNA sample preparations. Sequencing libraries were generated using NEBNext® UltraTM RNA Library Prep Kit for Illumina® (NEB, Ipswich, MA, USA) following the manufacturer's recommendations and index codes were added to attribute sequences to each sample. Briefly, mRNA was purified from total RNA using poly-T oligo-attached magnetic beads. Fragmentation was carried out using divalent cations under elevated temperature in NEBNext First Strand Synthesis Reaction Buffer (5X). First-strand cDNA was synthesized using random hexamer primer and M-MuLV Reverse Transcriptase (RNase H-). Second strand cDNA synthesis was subsequently performed using DNA Polymerase I and RNase H. Remaining overhangs were converted into blunt ends via exonuclease/polymerase activities. After adenylation of 3' ends of DNA fragments, NEBNext Adaptor with hairpin loop structure were ligated to prepare for hybridization. To select cDNA fragments of ~150-200 bp in length, the library fragments were purified with AMPure XP system (Beckman Coulter, Beverly, USA). Then 3 μ L USER Enzyme (NEB, Ipswich, MA, USA) was used with size-selected, adaptor-ligated cDNA at 37 °C for 15 minutes followed by five minutes at 95 °C before PCR. PCR was performed with Phusion High-Fidelity DNA polymerase, Universal PCR primers and Index (X) Primer. The resulting PCR products were purified (AMPure XP system) and library quality was assessed on the Agilent Bioanalyzer 2100 system. Sequencing was carried out on an Illumina NovaSeq® 6000 (NovoGene; Illumina, San Diego, California, United States), using paired-end 2 \times 150 bp chemistry at a depth of 20 million reads per sample.

Sequencing quality was examined using FastQC(12) (v0.11.2), and adaptors were subsequently trimmed with Trimmomatic(13) (v0.36) (in paired-end mode, parameters: TruSeq3-PE-2.fa:2:30:7:8:true LEADING:10 TRAILING:10 SLIDINGWINDOW:5:20 MINLEN:36). Trimmed fastq files were subsequently re-examined with FastQC again to ensure efficient adaptor and quality trimming. Reads were aligned to the hg38 genome (obtained through UCSC) with STAR(14) (v2.7.8a) using the default parameters only for quality control purposes. Aligned reads were subjected to the following quality control metrics: (i) duplication rate assessment through Picard MarkDuplicates (v.2.18), (ii) rRNA content, genomic read distribution and 3'-5' bias through RNA-SeQC (v.2.4.0). These metrics, together with the read alignment summary from STAR, were summarized using MultiQC(15) (v1.9), resulting in the removal of one sample on the basis of abnormal read duplication rates (n=1 mild disease). Gene quantification was performed using Salmon(16) (v1.3.0) with adaptor-trimmed reads and gene annotation obtained from GENCODE version 31 (parameters: --gcBias --validateMappings). Gene-level read counts were imported into R using R package Tximport(13)(v.1.18.0) with transcript length adjustment. Read counts were first normalized by DESeq2(17) (v.1.30.0), with the svaseq() function from R package sva (v.3.38.0)(18) being used to subsequently estimate hidden batch effects with n.sv=4. Sva-normalized read counts were obtained by regressing out covariates using the following function: <https://github.com/LieberInstitute/jaffelab/blob/master/R/cleaningY.R> and were subsequently used for PCA analysis and sample correlation analysis. R package fgsea (v 1.16.0)(19) was used to perform Gene Set Enrichment Analysis with genes ranked based on fold change between the compared conditions. Gene set files were obtained from the Molecular Signatures Database (MSigDB, v7.0) and only C2: curated gene sets and C5: ontology gene sets were used. Pathway enrichment using significantly differentially expressed genes were performed with R package gProfilerR (v.0.7.0)(20). R scripts and essential data files to reproduce the mRNA-seq analysis can be found in the Online Data Files.

HTG EdgeSeq MicroRNA (miRNA) Whole Transcriptome Assay (WTA) from Plasma: Lysis of t_0 -1 plasma aliquots was facilitated by combining 30 μ L plasma with equivalent (v/v) amounts of HTG Plasma Lysis Buffer (HTG Molecular, Tucson, AZ, USA) as well as $1/10^{\text{th}}$ (v/v) amounts of Proteinase K (HTG Molecular, Tucson, AZ, USA). The mixture was subsequently incubated for three hours at 50°C shaking at 1,400 rpm. From each prepared sample, 35 μ L were added per well to a 96-well sample plate. Human fetal brain RNA was added to one well at 25 ng/well to serve as an internal control. Samples were run on an HTG EdgeSeq Processor using the HTG EdgeSeq miRNA WTA (HTG Molecular, Tucson, AZ, USA) to facilitate nuclease protection, whereby a pre-selected miRNA population is protected with proprietary protection probes, followed by degradation of all non-hybridized probes and non-targeted RNA by S1 nuclease. Following the processing, samples were individually barcoded (using a 16-cycle PCR reaction), individually purified using AMPure XP beads (Beckman Coulter, Brea, CA, USA), and quantified using a KAPA Library Quantification kit (KAPA Biosystem, Wilmington, MA, USA). The library was sequenced on a NextSeq (Illumina, Inc., San Diego, CA) using a V3 150-cycle kit with two index reads. PhiX (Roche, Mississauga, ON, CAN) was spiked into the library at 5%; this spike-in control is standard for Illumina sequencing libraries. Data were returned from the sequencer in the form of demultiplexed FASTQ files, with one file per original well of the assay. The HTG EdgeSeq Parser (v. 5.0.535.3181, HTG Molecular, Tucson, AZ, USA) was used to align the FASTQ files to the probe list to collate the data. Data were provided as data tables of raw, quality control (QC) raw, counts per million, and median normalized.

HTG EdgeSeq MiRNA Analysis: Samples were initially analyzed using three QC metrics: (i) QC0, examining degraded sample; cut-off of positive % $\geq 14\%$ as failure, (ii) QC1, insufficient read depth; read depth $\leq 500\text{k}$ as failure, (iii) QC2, minimal expression variability; relative standard deviation of reads ≤ 0.08 as failure. Of the 156 samples sent for sequencing, 12 exhibited failure at the level of QC2 and were excluded from all downstream analyses. Normalization of miRNA expression data on the remaining samples was performed using DeSeq2(17) (v. 1.14.1) in the HTG reveal software (v.3.0.0, HTG Molecular, Tucson, AZ, USA). MiRNAs were considered detectable if they had expression levels of >5 counts per million in more than half of our samples. Expression counts were logarithmically scaled (\log_{10}) for data visualization.

Simultaneous Multiplexed *In Vitro* Immunofluorescence: Cells were fixed in ice-cold methanol (Sigma-Aldrich, St Louis, MO) for five minutes at room temperature. Blocking was subsequently conducted in a solution containing 1% BSA (w/v, BioShop, Burlington, ON, CAN), 22.52 mg/mL glycine (ThermoFisher, Waltham, MA, USA), and PBS^{-/-} with Tween 20 (PBST, PBS^{-/-} + 0.1% Tween 20) for 30 minutes. Immunostaining was then conducted with diluted antibody (Online Table II) in 1% BSA PBST for 16 hours at 4°C. The cells were subsequently washed three times in PBS^{-/-}, five minutes for each wash, and re-incubated with the concordant secondary antibodies in 1% BSA for one hour at room temperature in the dark. Post-washing, cells were counterstained with Vectashield Antifade Mounting Media with 4',6-diamidino-2-phenylindole ([DAPI], Vector Laboratories, Burlingame, CA, USA) and mounted with coverslips (VWR International, Mississauga, ON, CAN) and sealed with nail polish. Confocal images were taken using an Olympus Fluoview 1000 Confocal microscope Olympus IX81 inverted stand (Olympus, CA, USA). Fluorochromes were excited using the following wavelengths: 405 nm for DAPI, 473 nm for Alexa Fluor 488, and 559 nm for Alexa Fluor 568 (ThermoFisher, Waltham, MA, USA). A 20X/0.75NA UPLSAPO super apochromat objective was used to take the 20X images while a Plan Apo 40x/1.35 NA oil immersion objective was utilized for the 40X images. Image processing was done with the FV10-ASW 4.2 Viewer (Olympus, CA, USA). Image intensities were calculated using FIJI(21) (v2.1.0/1.53c), by examining the integrated density and dividing that by the number of DAPI positive cells within each image.

SUPPLEMENTAL TABLES

Online Table I. Key Reagents and Resources

REAGENT or RESOURCE	SOURCE	IDENTIFIER
Enzyme-Linked Immunosorbent Assays		
Ang-2, sE-SEL, sICAM-1, sVCAM-1	ProteinSimple	Cat#: SPCKC-PS-004112
ET-1, IL-6, IL-8, sTREM-1	ProteinSimple	Cat#: SPCKC-PS-004111
High Sensitivity Cardiac Troponin I	Biomatik	Cat#: EKV09460
MPO	R&D Systems	Cat#: DMYE00B
sSlit2	Elabscience	Cat#: E-EL-H0931
SARS-CoV-2 Spike (Trimer) Ig Total	ThermoFisher	Cat#: BMS2323
Chemicals, Peptides, and Recombinant Proteins		
Bovine Serum Albumin	BioShop	Cat#: ALBC0100
Dexamethasone	BioShop	Cat#: DEX002.100
EBM-2 TM	Lonza	Cat#: 00190860
EGM-2 TM Bullet Kit	Lonza	Cat#: CC-3162
FITC-40kDa Dextran	Sigma-Aldrich	Cat#: 53379
Gelatin	ThermoFisher	Cat#: S006100
Hanks' Balanced Salt Solution	Thermo	Cat#: 14025092
Methanol	Sigma-Aldrich	Cat#: 322415
Mounting Medium with DAPI	Vectashield	Cat#: H-1200
Nangibotide	LifeTein	Cat#: Custom Order
Phosphate Buffered Saline	Gibco	Cat#: LS10010023
Q-Peptide	Genscript	Cat#: SC1208
RIPA Buffer (10x)	EMD Millipore	Cat#: 20-188
Slit2-N	PreproTech	Cat#: 150-11
Thrombin	Sigma-Aldrich	Cat#: 10602400001
TNF α Human	Sigma-Aldrich	Cat#: T0157-10UG
TritonX-100	Sigma-Aldrich	Cat#: T8787
Tween20	ThermoFisher	Cat#: 003005
UltraPure Glycine	ThermoFisher	Cat#: 15527013
Commercial Components		
Vacutainer EDTA Tubes	BD Diagnostics	Cat#: 367525
E-16 Plates	Agilent	Cat#: 300601150
Black 96 well assay plate	Sigma-Aldrich	Cat#: M0312
Costar TC-Treated Plates (24-well)	Millipore Sigma	Cat#: CLS3527-100EA
Coverslips No. 1, 24x50mm	VWR	Cat#: 4839081
DNA LoBind Tubes	Eppendorf	Cat#: 22431021
Mycoplasma Contamination Kit	ThermoFisher	Cat#: 4460623
RNeasy Plus Micro Kit	Qiagen	Cat#: 74034
Transwell Filters (12-well plate)	Corning	Cat#: 3462
NEBNext® Ultra TM Library Prep Kit	New England Biolabs	Cat#: E7645S
AMPure XP Beads	Beckman Coulter	Cat#: A63880
Uracil-Specific Excision Reagent	New England Biolabs	Cat#: M5505S
Phusion High-Fidelity DNA polymerase	ThermoFisher	Cat#: F-530XL
HTG Lysis Buffer	HTG Diagnostic's	Cat#: SPP-Mi-04
KAPA Library Quantification Kit	Roche	Cat#: 07960140001
PhiX	Illumina	Cat#: FC-110-3001
Cells		
pHUEVC	Lonza	Cat#: C2519A Lot: 661173

Cat#: C2519A Lot: 636514		
Key Instruments		
Cytation 5	Biotek	Serial Number: 16041913
DS-11 + Spectrophotometer	DeNovix	Serial Number: 760419B
ELLA	ProteinSimple	Serial Number: ELLA-16080112
NextSeq	Illumina	Serial Number: A0877
NovoSeq6000	Illumina	Serial Number:
Fluoview 1000 Confocal IX81	Olympus	Serial Number: 8B03839
Real-Time Cell Analysis	xCelligence	Serial Number: 3211107167871
Software and Algorithms		
Adobe Illustrator	https://www.adobe.com/products/illustrator.html (Ver. 25.4.1)	
DeSeq2 (17)	https://bioconductor.org/packages/release/bioc/html/DESeq2.html (Ver. 3.12)	
EdgeSeq Parser	HTG Diagnostics (Ver. 5.0.525.3181)	
FastQC (12)	https://www.bioinformatics.babraham.ac.uk/projects/fastqc/ (Ver. 0.11.2)	
FIJI (21)	https://imagej.net/software/fiji/downloads (Ver. 2.1.0/1.53c)	
FV10-ASW Viewer	https://www.olympus-lifescience.com/es/support/downloads (Ver. 4.2b)	
Fgsea (v 1.16.0) (19)	https://github.com/ctlab/fgsea#~:text=fgsea%20is%20an%20R%2Dpackage,level%20split%20Monte%2DCarlo%20scheme. (Ver. 1.16.0)	
gProfilerR	https://biit.cs.ut.ee/gprofiler/ (Ver. 0.7.0)	
HTG Reveal	HTG Diagnostics (Ver. 3.0.0)	
MarkDuplicates	http://broadinstitute.github.io/picard (Ver. 2.18)	
miRPath (22)	http://snf-515788.vm.okeanos.grnet.gr/ (Ver. 3.0)	
MultiQC (15)	https://github.com/ewels/MultiQC (Ver. 1.9.0)	
Prism 9	GraphPad Software ((Ver. 9.0.0 (86))	
Python	https://www.python.org/ (Ver. 3.8.8)	
Qualimap (23)	https://github.com/EagleGenomics-cookbooks/QualiMap (Ver. 2.2.1)	
R (24)	https://www.r-project.org/ (Ver. 4.0.3)	
RNA-SeQC (25)	https://github.com/getzlab/rnaseqc (Ver. 2.4.0)	
RTCA	xCelligence (Ver. 2.0.0)	
Salmon (16)	https://combine-lab.github.io/salmon/ (Ver. 1.3.0)	
Scikit-learn (26)	https://scikit-learn.org/stable/ (Ver. 0.24.1)	
Simple Plex Runner	Protein Simple (Ver. 3.7.1.12)	
STAR (14)	https://github.com/alexdobin/STAR (Ver. 2.7.8)	
Sva (18)	https://bioconductor.org/packages/release/bioc/html/sva.html (Ver.3.38.0)	
Trimmomatic (13)	https://bioconductor.org/packages/release/bioc/html/tximport.html (Ver. 0.36)	
Tximport (13)	https://bioconductor.org/packages/release/bioc/html/tximport.html (Ver. 3.12)	
Deposited Data		
Plasma microRNA transcriptome	Human participants	GSE: 178246
Messenger RNA sequencing	Pooled HUVECs	GSE: 178331

Abbreviations: Ang-2 = Angiopoietin-2; Cat# = Catalogue number; DAPI = 4',6-diamidino-2-phenylindole; EBM = Endothelial basal media; EGM = Endothelial growth medium; ET-1 = Endothelin-1; FITC = Fluorescein isothiocyanate; pHUVEC = Pooled human umbilical vein endothelial cells; Ig = Immunoglobulin; IL = Interleukin; kDa = Kilodalton; MPO = Myeloperoxidase; RIPA = Radioimmunoprecipitation assay buffer; SARS-CoV-2 = Severe acute respiratory syndrome coronavirus; sICAM = Soluble intercellular adhesion molecule-1; sSLIT-2 = Soluble slit guidance ligand 2; sTREM-1 = Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1; TC = Tissue culture; TNF α = Tumor necrosis factor alpha.

Online Table II. Antibodies

ANTIGEN	HOST	SOURCE	IDENTIFER
Human Vascular Endothelial-Cadherin	Mouse	R&D Systems	Cat#: 9381; RRID:AB_2260374
Polyclonal Claudin-5	Rabbit	Invitrogen	Cat#: 34-1600; RRID:AB_2533157
Alexa Fluor 488 Anti-rabbit	Goat	ThermoFisher	Cat#: A11008; RRID:AB_143165
Alexa Fluor 568 Anti-mouse	Donkey	ThermoFisher	Cat#: A10042; RRID:AB_2534017

Abbreviations: Cat# = Catalogue number; Research Resource Identifier (Antibody) = RRID:AB.

Online Table III. Admission Clinical Laboratory Findings

Characteristics*	Disease Severity					P-value
	Mild Negative (n=30)	Mild (n=27)	Moderate (n=39)	Severe (n=76)	Severe Negative (n=69)	
Complete Blood Count						
White blood cells – million/mm ³						
Median (IQR)	8.0 (6.1-10.0)	7.0 (5.2-11.2)	7.0 (5.1-9.1)	10.8 (8.3-13.8)	11.9 (7.4-16.9)	<0.001
Distribution – no./total no. (%)						
<4	3/30 (10.0)	4/27 (14.8)	4/36 (11.1)	2/75 (2.7)	4/69 (5.8)	0.1954
>11	3/30 (10.0)	7/27 (25.9)	6/36 (16.7)	35/75 (46.7)	36/69 (52.2)	<0.0001
Neutrophils – per mm ³						
Median (IQR)	5.1 (3.8-7.0)	4.7 (3.6-8.2)	4.7 (3.2-6.5)	9.0 (6.0-12.0)	9.5 (6.0-15.0)	<0.001
<3,000	4/30 (13.3)	5/27 (18.5)	8/36 (22.2)	1/75 (1.3)	3/68 (4.4)	0.0009
>5,800	11/30 (36.7)	10/27 (37.0)	12/36 (33.3)	22/75 (29.3)	14/68 (20.6)	0.3613
Lymphocytes – per mm ³						
Median (IQR)	1.4 (1.0-1.9)	1.3 (0.9-2.1)	1.2 (0.8-1.7)	0.9 (0.6-1.2)	0.9 (0.6-1.5)	0.0023
Distribution – no./total no. (%)						
<1,500	16/30 (53.3)	17/27 (63.0)	25/36 (69.4)	61/75 (81.3)	52/68 (76.5)	0.0343
>3,000	1/30 (3.3)	1/27 (3.7)	2/36 (5.6)	3/75 (4.0)	4/68 (5.9)	0.9693
Monocytes – per mm ³						
Median (IQR)	0.6 (0.4-0.8)	0.6 (0.4-0.8)	0.6 (0.3-0.7)	0.5 (0.3-0.8)	0.8 (0.3-1.0)	0.5077
Distribution – no./total no. (%)						
<300	1/30 (3.3)	4/27 (14.8)	5/36 (13.9)	14/75 (18.7)	21/68 (30.9)	0.0180
>500	16/30 (53.3)	15/27 (55.6)	18/36 (50.0)	32/75 (42.7)	42/68 (61.8)	0.2459
Platelet count – thousand/mm ³						
Median (IQR)	232 (179-292)	247 (177-273)	246 (179-292)	227 (154-313)	217 (151-268)	0.6530
Distribution – no./total no. (%)						
<150,000	2/30 (6.7)	4/27 (14.8)	3/36 (8.3)	17/75 (22.7)	17/69 (24.6)	0.1117
>400,000	2/30 (6.7)	1/27 (3.7)	5/36 (13.9)	6/75 (8.0)	5/69 (7.2)	0.6457
Cardiac Laboratory Results						
hs-cTnI – ng/mL, median (IQR)	2.50 (0.92-6.20)	0.92 (0.92-6.50)	1.90 (0.92-7.66)	16.4 (7.42-56.0)	23.0 (10.0-62.5)	<0.0001
Coagulation Laboratory Results						
International normalized ratio						
<0.9	0/11 (0.0)	0/5 (0.0)	0/17 (0.0)	0/74 (0.0)	0/65 (0.0)	1.0000
>1.1	5/11 (45.5)	4/5 (80.0)	12/17 (70.6)	47/74 (63.5)	45/65 (69.2)	0.5329
Partial thromboplastin time						
<25	6/8 (75.0)	2/5 (40.0)	2/14 (14.3)	0/74 (0.0)	0/64 (0.0)	<0.0001
>40	0/8 (0.0)	1/5 (20.0)	3/14 (21.4)	10/74 (13.5)	4/64 (6.3)	0.2930
Supplemental Laboratory Results						
Albumin, g/liter						
<35	2/5 (40.0)	4/9 (44.4)	8/17 (47.1)	59/74 (79.7)	44/69 (63.8)	0.0136
>50	0/5 (0.0)	0/9 (0.0)	0/17 (0.0)	1/74 (1.4)	0/69 (0.0)	0.8513
Alanine aminotransferase, >40 U/liter	1/14 (7.1)	5/18 (27.8)	3/24 (12.5)	25/72 (34.7)	17/66 (25.8)	0.0791
Aspartate aminotransferase, >40 U/liter	3/14 (21.4)	5/18 (27.8)	4/24 (16.7)	41/72 (56.9)	26/66 (39.4)	0.0019
Creatine kinase, ≥200 U/liter	2/13 (15.4)	2/5 (40.0)	3/36 (8.3)	54/74 (73.0)	35/63 (55.6)	<0.0001
Creatinine, ≥133 μmol/liter	1/30 (3.3)	4/27 (14.8)	6/36 (16.7)	23/75 (30.7)	20/69 (29.0)	0.0149
D-dimer, ≥0.5 mg/liter	3/8 (37.5)	9/10 (90.0)	10/12 (83.3)	26/49 (53.1)	8/56 (14.3)	<0.0001
Lactate Dehydrogenase – μkatal/liter						
<2.34	4/8 (50.0)	12/15 (80.0)	8/12 (66.7)	4/38 (10.5)	1/60 (1.7)	<0.0001
>4.68	0/8 (0.0)	0/15 (0.0)	1/12 (8.3)	5/38 (13.2)	11/60 (18.3)	0.2629
Total bilirubin, >17.1 μmol/liter	2/13 (15.4)	1/16 (6.3)	0/17 (0.0)	17/71 (23.9)	12/68 (17.6)	0.1047
Minerals, median (IQR) – mmol/liter						
Potassium	3.0 (3.8-4.4)	3.9 (3.7-4.3)	4.1 (3.7-4.4)	4.2 (3.8-4.3)	4.1 (3.7-4.8)	0.4223
Sodium	139 (137-140)	137 (136-141)	138 (135-140)	138 (135-143)	137 (135-140)	0.5582

* Summary statistics, where n-value is not provided, are based on the full population indicated in the column heading.

Bolded log-ranked P values are significant (P values <0.05). Abbreviations: hs-cTnI = High sensitivity cardiac troponin I; IQR = Interquartile range. Percentages may not add up to 100% due to rounding.

Online Table IV. Association of Baseline Clinical Characteristics to Mortality Amongst All Admitted COVID-19 Patients.

Variables (Baseline)	Unadjusted Log-Rank P Value
Comorbidity - Gout	0.003
MiRNA - hsa-miR-6080	0.009
Comorbidity - Coronary artery disease	0.010
Protein - Ang-2	0.015
MiRNA - hsa-miR-199a-3p	0.017
MiRNA - hsa-miR-4793-5p	0.027
MiRNA - hsa-miR-181a-5p	0.028
MiRNA - hsa-miR-mir-30b-5p	0.035
MiRNA - hsa-miR-6750-5p	0.036
Clinical - Age of patient at hospital admission	0.046
Protein - MPO	0.079
History of stroke	0.115
Comorbidity - Heart failure	0.149
Comorbidity - Obesity	0.151
Protein - IL-8	0.156
Comorbidity - CKD	0.164
MiRNA - hsa-miR-301a-3p	0.167
Protein - sVCAM-1	0.200
MiRNA - hsa-miR-mir-30c-5p	0.200
CV Medication - Anticoagulant	0.230
Clinical - Male	0.240
History of dyslipidemia	0.250
Comorbidity - OSA	0.250
MiRNA - hsa-miR-146a-5p	0.260
MiRNA - hsa-miR-miR-30c-5p	0.260
Protein - IL-6	0.280
MiRNA - hsa-miR-1	0.280
Comorbidity - Malignancy	0.280
History of cardiac procedure/surgery	0.300
Comorbidity - GERD	0.310
History of any prior CV procedure	0.310
History of arrhythmia	0.320
Protein - sTREM-1	0.340
Protein - E-Selectin	0.350
Comorbidity - Renal disease	0.360
MiRNA - hsa-miR-mir-339-3p	0.380
Comorbidity - Immunocompromised	0.390
Protein - sICAM-1	0.400
CV medication - Number of medications	0.430
CV medication - ARB	0.430
MiRNA - hsa-miR-mir-4706	0.440
Other CV condition - PFO	0.440
Comorbidity - Hypertension	0.480
History of myocardial infarction	0.490
Other CV condition - Aortic aneurysm	0.490
CV medication - Statin	0.510
Comorbidity - Diabetes	0.520
CV medication - ACE inhibitor	0.520
CV medication - CCB	0.530
Other CV condition - TIA	0.550
Other CV condition - Ischemic heart disease	0.560
Comorbidity - Vascular disease	0.620
Comorbidity - Pneumonia	0.660
Other CV condition - PAD	0.680
Comorbidity - COPD	0.690
MiRNA - hsa-miR-mir-26a-5p	0.720
Clinical - Patient ethnicity	0.730
CV medication - Beta blocker	0.740
MiRNA - hsa-miR-mir-30d-5p	0.750
CV medication - Diuretics	0.760
Comorbidity - CF	0.800
Other CV condition - Endocarditis	0.800
Other CV condition - Dilatated aortic root	0.800
Other CV condition - Becker's muscle dystrophy	0.800
Comorbidity - Asthma	0.850

Comorbidity - Valvular heart disease	0.930
Other CV condition - VTE	0.950
Other CV condition - NSTEMI	1.000

* Gout (n=3) was a small number of observations.

Bolded log-ranked P values are significant (P values <0.05). Abbreviations: ACE = Angiotensin-converting-enzyme inhibitors; Ang-2 = Angiopoietin-2; ARB = Angiotensin II receptor blockers; CCB = Calcium channel blocker; CKD = Chronic Kidney Disease; CF = Cystic Fibrosis; COPD = Chronic obstructive pulmonary disease; GERD = Gastroesophageal reflux disease; IL = Interleukin; miR = MicroRNA; MPO = Myeloperoxidase; NSTEMI = Non-ST-elevation myocardial infarction; OSA = Obstructive Sleep Apnea; PAD = Peripheral arterial disease; PFO = Patent foramen ovale; sICAM, Soluble intercellular adhesion molecule-1; sTREM-1, Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1; TIA = Transient ischemic attack; VTE = Venous thromboembolism.

Online Table V. Intensive Care Unit-Level Clinical Characteristics and Clinical Outcomes

Characteristics*	Severe (n=76)	Severe Negative (n=69)	P-value
Demographics			
Age, median (IQR) – yr.	61 (52-71)	61 (51-72)	0.5137
Distribution – no. (%)			
18-40 yr.	2 (2.6)	11 (15.9)	0.0071
41-64 yr.	45 (59.2)	26 (37.7)	0.0126
≥65 yr.	29 (38.2)	32 (46.4)	0.3999
Male Sex, no./total no. (%)	52/76 (68.4)	46/69 (66.7)	0.8602
BMI, median (IQR) [†]	28.1 (24.2-31.8)	25.7 (21.8-32.5)	0.3641
Obesity [‡] – no. (%)	38/61 (62.3)	35/69 (50.7)	0.2168
Length of hospital stay, median (IQR) – days [†]	13 (7-35)	8 (3-15)	<0.0001
Max temperature, median (IQR), °C	36.8 (36.3-37.8)	36.4 (35.8-37.1)	0.0797
Mix temperature, median (IQR), °C	36.7 (36.3-37.3)	36.3 (35.6-37.1)	0.0815
Illness Severity			
APACHE II	21 (16-27)	20 (16-29)	0.8521
SOFA score			
Median, (IQR)	9 (4-10)	7 (3-10)	0.8112
≥2 – no. (%)	12/38 (31.6)	18/69 (26.1)	0.3674
≥6 – no. (%)	24/38 (63.2)	44/69 (63.8)	>0.9999
Respiratory metrics			
Proned – no. (%)	8/38 (21.1)	3/69 (4.3)	0.0158
Intubation – no. (%)	53/76 (69.7)	47/69 (68.1)	0.8565
NIV – mean, days [†]	3.12	1.5	0.2177
MIV – mean, days [†]	10.84	7.61	0.2215
ECMO – no. (%)	21/76 (27.6)	0 (0.0)	<0.0001
FiO ₂ , median % (IQR)	0.53 (0.50-0.63)	0.40 (0.30-0.50)	<0.0001
PF ratio, median (IQR), mm Hg/%	130 (103-177)	188 (146-289)	0.0003
PF ratio <300 mm Hg/%	26/26 (100.0)	38/49 (77.6)	0.0126
RR, no. (%), ≥22 breaths/min	22/38 (57.9)	33/65 (50.8)	0.5425
Cardiovascular metrics			
Heart rate, median (IQR)	80 (67-101)	82 (67-97)	0.6911
Blood Pressure, mmHg			
Mean (S.D.)	118.6 (19.26)	109.9 (19.25)	0.0334
Max Systolic, median (IQR),	165 (141-179)	140 (125-165)	0.0242
Max Diastolic, median (IQR),	75 (65-84)	69 (63-80)	0.1707
Therapies – no. (%)			
Intravenous antibiotics	11/38 (28.9)	20/69 (29.0)	>0.9999
Systemic glucocorticoids	13/38 (34.2)	37/67 (55.2)	0.2016
Remdesivir	1/38 (2.6)	0/69 (0.0)	0.3679
Fludrocortisone	0/38 (0.0)	5/69 (7.2)	0.1583
Outcomes – no. (%)			
Any secondary CV event	13/76 (17.1)	1/69 (1.4)	0.0013
ARDS during hospitalization	43/76 (56.6)	25/69 (36.2)	0.0195
Arrhythmia during hospitalization	8/76 (10.5)	0/69 (0.0)	0.0068

* Summary statistics, where n-value is not provided, are based on the full population indicated in the column heading.

[†] Body mass index is the weight in kilograms divided by the square of the height in meters.

[‡] Obesity is classified according to World Health Organization guidelines (i.e., >25).

Bolded log-ranked P values are significant (P values <0.05). Abbreviations: APACHE = Acute physiologic assessment and chronic health evaluation; ARDS = Acute respiratory distress syndrome; BMI = Body mass index; CV = Cardiovascular; ECMO = Extracorporeal membrane oxygenation; FiO₂ = Fraction of inspired oxygen; IQR = Interquartile range; MIV = Mechanical ventilation; NIV = Non-invasive ventilation; PF = Ratio of arterial oxygen partial pressure to fractional inspired oxygen; RR = Respiratory rate; S.D. = Standard deviation; SOFA = Sequential Organ Failure Assessment. Percentages may not add up to 100% due to rounding.


Online Table VI. Association of Baseline Clinical Characteristics to Mortality Amongst All Severe COVID-19 Patients.

Variables (Baseline)	Unadjusted Log-Rank P Value
Comorbidity – Gout	<0.001
MiRNA - miR-30b-5p	0.005
MiRNA - hsa-miR-199a-3p	0.006
MiRNA - hsa-miR-181a-5p	0.010
Protein - VCAM-1	0.012
MiRNA - hsa-miR-339-3p	0.017
MiRNA - hsa-miR-mir-30e-5p	0.017
Protein - Ang-2	0.020
MiRNA - hsa-miR-146a-5p	0.030
MiRNA - hsa-miR-30c-5p	0.035
Comorbidity - CKD	0.038
History of coronary artery disease at baseline	0.041
MiRNA - hsa-miR-6080	0.042
MiRNA - hsa-miR-6750-5p	0.051
MiRNA - hsa-miR-4793-5p	0.060
Protein - ICAM-1	0.096
MiRNA - hsa-miR-301a-3p	0.103
CV medication - Anticoagulant	0.103
Clinical - Age of patient at hospital admission	0.122
Comorbidity - OSA	0.126
Clinical - Male	0.130
MiRNA - hsa-miR-4706	0.131
Comorbidity - Stroke	0.184
Protein - IL-8	0.220
Comorbidity - Heart failure	0.220
MiRNA - hsa-miR-26a-5p	0.230
Protein - sTREM-1	0.230
History of arrhythmia	0.240
Comorbidity - Obesity	0.260
Comorbidity - COPD	0.260
MiRNA - hsa-miR-1	0.270
History of any prior CV procedure	0.310
Protein - IL-6	0.320
Protein - MPO	0.320
Comorbidity - GERD	0.320
Comorbidity - Dyslipidemia	0.380
Other CV condition - TIA	0.380
Comorbidity - Peripheral vascular disease	0.410
Comorbidity - Immunocompromised	0.410
Comorbidity - Valvular heart disease	0.420
Other CV condition - PFO	0.440
MiRNA - hsa-miR-30d-5p	0.460
CV medication - Statin	0.480
CV medication - Diuretics	0.540
CV medication - Number of medications	0.550
CV medication - ARB	0.550
Comorbidity - Renal disease	0.560
Other CV condition – Ischemic heart disease	0.560
Comorbidity - Hypertension	0.590
CV medication - CCB	0.590
History of myocardial infarction	0.600
Comorbidity - Malignancy	0.600
Other CV condition - Aortic aneurysm	0.620
CV medication - ACE inhibitor	0.650
Comorbidity - Pneumonia	0.660
Clinical - Patient ethnicity	0.730
Comorbidity - Diabetes	0.780
Protein - sE-Selectin	0.810
Comorbidity - Asthma	0.810
Other CV condition - VTE	0.840
CV medication - Beta blocker	0.880


* Gout (n=3) was a small number of observations.

Bolded log-ranked P values are significant (P values <0.05). Abbreviations: ACE = Angiotensin-converting-enzyme inhibitors; Ang-2 = Angiopoietin-2; ARB = Angiotensin II receptor blockers; CCB = Calcium channel blocker; CKD = Chronic Kidney Disease; CF = Cystic Fibrosis; COPD = Chronic obstructive pulmonary disease; GERD = Gastroesophageal reflux disease; IL = Interleukin; miR = MicroRNA; MPO = Myeloperoxidase; OSA = Obstructive Sleep Apnea; PFO = Patent foramen ovale; sICAM, Soluble intercellular adhesion molecule-1; sTREM-1, Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1; TIA = Transient ischemic attack; VTE = Venous thromboembolism.


SUPPLEMENTAL FIGURES



Online Figure I; Related to Methods and Table 1: Flow diagram of patients enrolled between the COLOBILI Study (St. Michael's Hospital) and the COVID Study (University Health Network). Abbreviations: COVID-19 = Coronavirus disease 2019; ED = Emergency department; ICU = Intensive care unit; RTI = Respiratory tract infection; SMH = St. Michael's Hospital; UHN = University Health Network.




Online Figure II; Related to Methods and Table 1: Spike (trimer) antigen serology testing from patients having a negative SARS-CoV-2 polymerase chain reaction result. The data were analyzed using built-in low and high thresholds, whereby the area between the indeterminate-low and indeterminate-high constitutes an ambiguous result; n=80. Positive control is indicated on the graph (red circle). The graph depicts averaged values of independent technical duplicates data with center bars representing the mean and error bars representing standard deviation (\pm S.D.).




Online Figure III; Related to Figure 1: The association of coronary artery disease with mortality characterized in terms of proportion of deceased patients stratified by status. The data were analyzed using log-rank testing. Abbreviations: CAD = Coronary artery disease.

Entire Cohort (SARS-CoV-2 Negative and Positive)




Online Figure IV; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers amongst the entire cohort (SARS-CoV-2 negative and positive populations). Abbreviations: Ang-2 = Angiopoietin-2; ET-1 = Endothelin-1; Hs-CTnI = High-sensitivity cardiac troponin I; IL = Interleukin; MPO = Myeloperoxidase; sICAM = Soluble intercellular adhesion molecule-1; sTREM-1 = Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1; WBC = White blood cells.

COVID-19 – Mild




Online Figure V; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the mild COVID-19 subgroup. Abbreviations: Ang-2 = Angiopoietin-2; ET-1 = Endothelin-1; Hs-cTnI = High-sensitivity cardiac troponin I; IL = Interleukin; MPO = Myeloperoxidase; sICAM = Soluble intercellular adhesion molecule-1; sTREM-1 = Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1; WBC = White blood cells.

COVID-19 – Moderate




Online Figure VI; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the moderate COVID-19 subgroup. Abbreviations: Ang-2 = Angiopoietin-2; ET-1 = Endothelin-1; Hs-cTnI = High-sensitivity cardiac troponin I; IL = Interleukin; MPO = Myeloperoxidase; sICAM = Soluble intercellular adhesion molecule-1; sTREM-1 = Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1; WBC = White blood cells.

COVID-19 – Severe




Online Figure VII; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the severe COVID-19 subgroup. Abbreviations: Ang-2 = Angiopoietin-2; ET-1 = Endothelin-1; Hs-cTnI = High-sensitivity cardiac troponin I; IL = Interleukin; MPO = Myeloperoxidase; sICAM = Soluble intercellular adhesion molecule-1; sTREM-1 = Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1; WBC = White blood cells.

SARS-CoV-2 Negative – Mild




Online Figure VIII; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the mild SARS-CoV-2 negative subgroup. Abbreviations: Ang-2 = Angiopoietin-2; ET-1 = Endothelin-1; Hs-cTnI = High-sensitivity cardiac troponin I; IL = Interleukin; MPO = Myeloperoxidase; sICAM = Soluble intercellular adhesion molecule-1; sTREM-1 = Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1; WBC = White blood cells.

SARS-CoV-2 Negative – Severe




Online Figure IX; Related to Figure 2: Spearman correlations between t_{0-1} concentrations of biomarkers within the severe SARS-CoV-2 negative subgroup. Abbreviations: Ang-2 = Angiopoietin-2; ET-1 = Endothelin-1; Hs-cTnI = High-sensitivity cardiac troponin I; IL = Interleukin; MPO = Myeloperoxidase; sICAM = Soluble intercellular adhesion molecule-1; sTREM-1 = Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1; WBC = White blood cells.




Online Figure X; Related to Figure 2. Plasma Concentration of Endothelial Dysfunction and Inflammatory Markers at t_{0-1} . (a) ET-1, (b) sICAM-1, (c) sE-Selectin, (d) sTREM-1, (e) IL-6, (f) IL-8, (g) MPO, and (h) hs-cTnI stratified


among disease severity. Data shown are for all patients with an available t_{0-1} sample ($n=210$), with values representing the mean and error bars are (\pm S.D.). P values for multiple group comparisons were determined by Kruskal-Wallis test with Dunn's multiple comparisons test. Severe negative comparisons are only shown in reference to the concordant severe group; testing was conducted with all groups. The graph depicts averaged values of independent technical triplicates data points with center bars representing the mean and error bars representing standard deviation (\pm S.D.). Abbreviations: ET-1 = Endothelin-1; hs-cTnI = High-sensitivity cardiac troponin I; sICAM = Soluble intercellular adhesion molecule-1; IL = Interleukin; MPO = Myeloperoxidase; sTREM-1 = Soluble triggering receptor expressed on myeloid cells 1; N.S. = non-significant.



Online Figure XI; Related to Figure 2: Plasma Concentration of Endothelial Dysfunction and Inflammatory Markers at t_{0-1} and ability to discriminate survival in ICU patients. Severe COVID-19 patients and severe negative patients (i.e., SARS-CoV-2 negative) (a) Levels of Angiopoietin-2, (b) Endothelin-1, (c) sICAM-1, (d) sVCAM-1, (e) sE-Selectin, (f) sTREM-1, (g) IL-6, (h) IL-8, (i) MPO, and (j) hs-cTnI stratified among disease severity. Data shown are for all severe patients with an available t_{0-1} sample ($n=114$), with the center bars representing the mean and error bars representing standard deviation (\pm S.D.). P values for multiple group comparisons were determined by Kruskal-Wallis test with Dunn's multiple comparisons test. Abbreviations: Ang-2 = Angiopoietin-2; sICAM = Soluble intercellular adhesion molecule 1; IL = Interleukin; MPO = Myeloperoxidase; sTREM-1 = Soluble triggering receptor expressed on myeloid cells 1; sVCAM-1 = Soluble vascular cell adhesion molecule 1.




Online Figure XII; Related to Figure 2: Plasma Concentration of Endothelial Dysfunction and Immunological Markers at t_{0-1} in ICU patients. (a) Levels of Angiopoietin-2, (b) Endothelin-1, (c) sICAM-1, (d) sVCAM-1, (e) sE-Selectin, (f) sTREM-1, (g) IL-6, (h) IL-8, (i) MPO, and (j) hs-cTnI stratified among disease severity. Data shown are for all patients with an available t_{0-1} ($n=114$), t_{2-3} ($n=84$), and t_{4-5} ($n=44$), with center bars representing the mean and error bars representing standard deviation (\pm S.D.). P values for multiple group comparisons were determined by 2-way ANOVA with Sidak's multiple comparisons test. Abbreviations: sICAM = Soluble intercellular adhesion molecule 1; IL = Interleukin; MPO = Myeloperoxidase; sTREM-1 = Soluble triggering receptor expressed on myeloid cells 1; sVCAM-1 = Soluble vascular cell adhesion molecule 1.




Online Figure XIII; Related to Figure 2: Plasma Concentration of (a) IL-6 and (b) MPO, longitudinally between severe COVID-19 patients and severe SARS-CoV-2 negative patients. The error bars represent the mean and its 95% confidence intervals estimated using generalized estimating equation with an independent working correlation matrix. The standard errors were estimated using robust sandwich estimator. Abbreviations: IL = Interleukin; MPO = Myeloperoxidase.

A




B

C



D

E




F

Online Figure XIV; Related to Figure 3 and 4: Plasma MicroRNA Transcriptome Across the Disease Severity Subgroups. Volcano plots of differentially expressed miRNA between patient groups (a, c, e) with predicted KEGG terms (with enrichment score below and number of genes to the right) for pathways of deregulated microRNAs shown beside each corresponding region of the volcano plot (b, d, f). Data is displayed as FDR adjusted P values (Q values) vs the \log_2 fold change, with dashed lines are drawn to define restriction boundaries. See Supplementary Data Files III and IV for a full list of differentially expressed miRNA along with a full list of predicted KEGG pathways. Abbreviations: COVID-19 = Coronavirus Disease 2019; FDR = False discovery rate; KEGG = Kyoto Encyclopedia of Genes and Genomes; MiRNAs = MicroRNAs; SARS-CoV-2 = Severe acute respiratory syndrome coronavirus 2.


A

Moderate COVID-19 vs Mild SARS-CoV-2 Negative




B

Concordantly Regulated Pathways




C

Severe COVID-19 vs Mild SARS-CoV-2 Negative




D

Concordantly Regulated Pathways




Concordantly Regulated Pathways




E

Severe Negative vs Mild SARS-CoV-2 Negative




F

Concordantly Regulated Pathways



Concordantly Regulated Pathways



Online Figure XV; Related to Figure 3 and 4: Plasma MicroRNA Transcriptome Across the Disease Severity Subgroups. Volcano plots of differentially expressed miRNA between patient groups (a, c, e) with predicted KEGG terms (with enrichment score below and number of genes to the right) for pathways of deregulated microRNAs shown beside each corresponding region of the volcano plot (b, d, f). Data is displayed as FDR adjusted P values (Q values) vs the \log_2 fold change, with dashed lines are drawn to define restriction boundaries. See Supplementary Data Files III and IV for a full list of differentially expressed miRNA along with a full list of predicted KEGG pathways. Abbreviations: COVID-19 = Coronavirus Disease 2019; FDR = False discovery rate; KEGG = Kyoto Encyclopedia of Genes and Genomes; MiRNAs = MicroRNAs; SARS-CoV-2 = Severe acute respiratory syndrome coronavirus 2.






Figure XVI; Related to Figure 4: Feature importance of a machine learning model incorporating clinical data. All clinical metrics are at time of admission with preexisting conditions defined according to those listed in the methods. Abbreviations: BMI = Body mass index; CAD = Coronary Artery Disease; COPD = Chronic obstructive pulmonary disease.




Online Figure XVII; Related to Figure 4: Feature importance of a machine learning model incorporating both clinical data and protein expression metrics. All clinical metrics are at time of admission with preexisting conditions defined according to those listed in the methods. Abbreviations: Ang-2 = Angiopoietin-2; BMI = Body mass index; CAD = Coronary Artery Disease; ET-1 = Endothelin-1; Hs-cTnI = High-sensitivity cardiac troponin I; IL = Interleukin; MPO = Myeloperoxidase; sICAM = Soluble intercellular adhesion molecule-1; sTREM-1 = Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1 = Soluble vascular cell adhesion molecule-1.




Online Figure XVIII; Related to Figure 4: Association of Biomarkers with In-Hospital Mortality for Severe COVID-19 Patients. Univariable log hazard ratios of candidate microRNAs (a) hsa-miR-181a-5p, (b) hsa-miR-199a-3p, and c) hsa-miR-399-3p.


A




B



Online Figure XIX; Related to Figure 5: T₀₋₁ COVID-19 Patient Plasma Selectively Induces Acute Increases in Endothelial Permeability. (a) Permeability of pHUVEC monolayers was measured by 40 kDa FITC extravasation from the apical to the basolateral surface one-hour post-co-incubation. Treatment groups were normalized to the negative control. Center line represents the mean and error bars are (\pm S.D.). P values determined by one-way ANOVA test with Tukey's multiple comparisons test. Moderate vs negative: $**P=9.5 \times 10^{-3}$; Severe vs negative $****P=5.6 \times 10^{-8}$; Severe vs moderate (P): $***P=6.5 \times 10^{-4}$; Severe vs severe negative: $****P=9.8 \times 10^{-5}$. (b) Permeability of pHUVEC monolayers was measured by 40 kDa FITC extravasation from the apical to the basolateral surface six hours post-co-incubation. Treatment groups were normalized to the negative control. Center bars represent the mean and error bars represent the standard deviation (\pm S.D.). P were values determined by one-way ANOVA test with Tukey's multiple comparisons test. Severe vs negative: $*P=2.1 \times 10^{-2}$; Severe vs severe negative: $*P=3.2 \times 10^{-2}$; n=7-8 per group. Thrombin treatment was included as a barrier disrupting positive control.



Online Figure XX; Related to Figure 5: Correlation of t₀₋₁ Plasma Cardiovascular Biomarkers in COVID-19 positive patients to Induction of Endothelial Permeability. Pearson correlations between (a) hemolysis, (b) Ang-2, (c) hs-cTnI, (d) sE-Selectin, (e) ET-1, (f) sICAM-1, (g) IL-6, (h) IL-8, (i) sTREM-1, (j) sVCAM-1, and (k) MPO to the change in pHUVEC TEER after six-hours co-incubation; n=111 per correlation. Leak is defined through negative values on the x-axis. Abbreviations: Ang-2 = Angiopoietin-2; ET-1 = Endothelin-1; Hs-cTnI = High-sensitivity cardiac troponin I; IL = Interleukin; MPO = Myeloperoxidase; sICAM = Soluble intercellular adhesion molecule-1; sTREM-1 = Soluble triggering receptor expressed on myeloid cells-1; sVCAM-1= Soluble vascular cell adhesion molecule-1.



Online Figure XXI; Related to Figure 6: Endogenous sSlit2 is upregulated in severe COVID-19 patient plasma. (a) Endogenous sSlit2 at t_{0-1} across the severity of COVID-19 ($n=27-40$, severe vs negative, $*P=0.0279$). (b) Endogenous sSlit2 at longitudinal intervals in patients with severe COVID-19 ($n=14-38$). Center bars represent the mean and error bars represent the standard deviation (\pm S.D.). P values were determined by one-way ANOVA test with Tukey's multiple comparisons test.

SUPPLEMENTAL DATA FILE ANNOTATIONS

Supplementary Data File I. Quality control table for all RNA-sequencing experiments used in this study.

Supplementary Data File II. R documentation file for the analysis of the RNA-sequencing experiments.

Supplementary Data File III. R code for the analysis of the RNA-sequencing experiments.

Supplementary Data File IV. Full list of differentially expressed miRNA with pairwise comparisons between COVID-19 cohorts and the negative controls.

Supplementary Data File V. Full list of pathway enrichments for miRNA-sequencing experiment between COVID-19 cohorts and the negative controls.

Supplementary Data File VI. Full list of differentially expressed genes with pairwise comparisons between COVID-19 cohorts and the negative controls.

Supplementary Data File VII. Full list of pathway enrichments for mRNA-sequencing experiment between COVID-19 cohorts and the negative controls.

Supplementary Data File VIII. Gene set enrichment analysis for mRNA-sequencing experiment between COVID-19 cohorts and the negative controls.

SUPPLEMENTAL REFERENCES

1. Health. NIo. COVID-19 Treatment Guidelines Panel. Coronavirus Disease 2019 (COVID-19) Treatment Guidelines. 2021 [Available from: <https://www.covid19treatmentguidelines.nih.gov/>. 2022. March 5
2. Ontario PH. COVID-19 Provincial Testing Guidance Update 2020 [updated September 24, 2020. Available from: http://www.health.gov.on.ca/en/pro/programs/publichealth/coronavirus/docs/2019_testing_guidance.pdf. 2022. March 5
3. Ontario PH. Ontario COVID-19 Data Tool 2021 [Available from: <https://www.publichealthontario.ca/en/data-and-analysis/infectious-disease/covid-19-data-surveillance/covid-19-data-tool>. 2022. March 5
4. Organization WH. [Available from: <https://www.who.int/news-room/fact-sheets/detail/obesity-and-overweight>. 2022. March 5
5. Larvin M, McMahon M. APACHE-II score for assessment and monitoring of acute pancreatitis. *Lancet*. 1989;**334**(8656):201-5.
6. Ferreira FL, Bota DP, Bross A, Mélot C, Vincent J-L. Serial evaluation of the SOFA score to predict outcome in critically ill patients. *J Am Med Assoc*. 2001;**286**(14):1754-8.
7. Miklas JW, Dallabrida SM, Reis LA, Ismail N, Rupnick M, Radisic M. QHREDGS Enhances Tube Formation, Metabolism and Survival of Endothelial Cells in Collagen-Chitosan Hydrogels. *Plos One*. 2013;**8**(8):e72956.
8. Sherchan P, Huang L, Akyol O, Reis C, Tang J, Zhang JH. Recombinant Slit2 Reduces Surgical Brain Injury Induced Blood Brain Barrier Disruption via Robo4 Dependent Rac1 Activation in a Rodent Model. *Sci Rep*. 2017;**7**(1):746.
9. François B, Wittebole X, Ferrer R, Mira J-P, Dugernier T, Gibot S, et al. Nangibotide in patients with septic shock: a Phase 2a randomized controlled clinical trial. *Intensive Care Med*. 2020;**46**(7):1425-37.
10. Murata T, Hori M, Lee S, Nakamura A, Kohama K, Karaki H, et al. Vascular endothelium has a local anti-adenovirus vector system and glucocorticoid optimizes its gene transduction. *Arterioscler Thromb Vasc Biol*. 2005;**25**(9):1796-803.
11. DiStefano PV, Kuebel JM, Sarelius IH, Glading AJ. KRIT1 Protein Depletion Modifies Endothelial Cell Behavior via Increased Vascular Endothelial Growth Factor (VEGF) Signaling. *J Biol Chem*. 2014;**289**(47):33054-65.
12. Andrews S. FastQC: a quality control tool for high throughput sequence data: Babraham Bioinformatics, Babraham Institute, Cambridge, United Kingdom; 2010 [Available from: <https://www.bioinformatics.babraham.ac.uk/projects/fastqc/>. 2022. March 5
13. Sonesson C, Love MI, Robinson MD. Differential analyses for RNA-seq: transcript-level estimates improve gene-level inferences. *F1000Research*. 2015;**4**.
14. Dobin A, Davis CA, Schlesinger F, Drenkow J, Zaleski C, Jha S, et al. STAR: ultrafast universal RNA-seq aligner. *Bioinformatics*. 2013;**29**(1):15-21.
15. Ewels P, Magnusson M, Lundin S, Käller M. MultiQC: summarize analysis results for multiple tools and samples in a single report. *Bioinformatics*. 2016;**32**(19):3047-8.
16. Patro R, Duggal G, Love MI, Irizarry RA, Kingsford C. Salmon provides fast and bias-aware quantification of transcript expression. *Nat Methods*. 2017;**14**(4):417-9.
17. Love MI, Huber W, Anders S. Moderated estimation of fold change and dispersion for RNA-seq data with DESeq2. *Genome Biol*. 2014;**15**(12):1-21.

18. Leek JT, Storey JD. Capturing heterogeneity in gene expression studies by surrogate variable analysis. *PLoS Genet.* 2007;**3**(9):1724-35.
19. Korotkevich G, Sukhov V, Budin N, Shpak B, Artyomov M, Sergushichev A. Fast gene set enrichment analysis. bioRxiv; 2016.
20. Reimand J, Arak T, Adler P, Kolberg L, Reisberg S, Peterson H, et al. g:Profiler-a web server for functional interpretation of gene lists (2016 update). *Nucleic Acids Res.* 2016;**44**(W1):W83-9.
21. Schindelin J, Arganda-Carreras I, Frise E, Kaynig V, Longair M, Pietzsch T, et al. Fiji: an open-source platform for biological-image analysis. *Nat Methods.* 2012;**9**(7):676-82.
22. Vlachos IS, Zagganas K, Paraskevopoulou MD, Georgakilas G, Karagkouni D, Vergoulis T, et al. DIANA-miRPath v3.0: deciphering microRNA function with experimental support. *Nucleic Acids Res.* 2015;**43**(W1):W460-W6.
23. Okonechnikov K, Conesa A, García-Alcalde F. Qualimap 2: advanced multi-sample quality control for high-throughput sequencing data. *Bioinformatics.* 2016;**32**(2):292-4.
24. Team RC. R: A language and environment for statistical computing. 2020.
25. DeLuca DS, Levin JZ, Sivachenko A, Fennell T, Nazaire M-D, Williams C, et al. RNA-SeQC: RNA-seq metrics for quality control and process optimization. *Bioinformatics.* 2012;**28**(11):1530-2.
26. Pedregosa F, Varoquaux G, Gramfort A, Michel V, Thirion B, Grisel O, et al. Scikit-learn: Machine learning in Python. *Journal of Machine Learning research.* 2011;**12**:2825-30.